# TaKaRa Ex Premier™ DNA Polymerase Dye plus

**Code No. RR371S      Size:      1 ml**  
**(for 40 PCR reactions)**

## Description:

TaKaRa Ex Premier DNA Polymerase Dye plus is a novel PCR enzyme based on DNA polymerase derived from *Thermococcus* sp. This product improves the PCR success rate for GC or AT rich sequences, long targets, and crude samples that are difficult to amplify with general PCR reagents due to the optimization of enzymes and buffers. Furthermore, it has high-fidelity ability of  $\alpha$ -type DNA polymerases and can be used for cloning.

This product is 2X premix type that contains all components (enzyme, buffer, dNTP Mixture, etc.) and does not freeze at -20°C, so it is possible to perform PCR by adding primers and templates without thawing. In addition, this product also contains a gel loading dye (green) that is convenient for electrophoresis, therefore the PCR products amplified with this master mix can be loaded directly onto an agarose gel electrophoresis.

This enzyme is designed for hot start PCR using the monoclonal antibody, thus it is possible to prevent non-specific amplification due to mis-priming and/or formation of primer dimer before PCR reaction.

**Storage:** -20°C

Note: Gently mix well before use and centrifuge briefly.

## Composition of the PCR reaction mixture (total 50 $\mu$ l):

|  |  |
|--|--|
| TaKaRa Ex Premier DNA Polymerase Dye plus (2X) | 25 $\mu$ l (final conc. 1X)                  |
| primer 1                                       | 10 - 15 pmol (final conc. 0.2 - 0.3 $\mu$ M) |
| primer 2                                       | 10 - 15 pmol (final conc. 0.2 - 0.3 $\mu$ M) |
| Template                                       | See "Template" chapter                       |
| Sterile purified water                         | up to 50 $\mu$ l                             |

## PCR conditions:

|           |                         |             |
|-----------|-------------------------|-------------|
| 94°C      | 1 min <sup>*1</sup>     |             |
| ↓         |                         |             |
| 98°C      | 10 sec                  | ] 30 cycles |
| 55 - 60°C | 15 sec                  |             |
| 68°C      | 30 sec/kb <sup>*2</sup> |             |

[ For GC-rich targets or amplification of products  $\geq$  10 kb in length ]

|      |                         |             |
|------|-------------------------|-------------|
| 94°C | 1 min <sup>*1</sup>     |             |
| ↓    |                         |             |
| 98°C | 10 sec                  | ] 30 cycles |
| 68°C | 30 sec/kb <sup>*2</sup> |             |

\*1 When amplifying a GC-rich region or a long target, perform initial denaturation at 94°C for 1 min.

\*2 When amplifying from crude sample, set the extension time to 1 min/kb.

Note 1: Do not use inosine-containing primers with TaKaRa Ex Premier DNA Polymerase Dye plus.

Note 2: Please set the annealing temperatures according to the T<sub>m</sub> value of the primers.

## Template:

Recommended quantities of template DNA (for 50  $\mu$ l reaction):

|                            |                 |                         |
|----------------------------|-----------------|-------------------------|
|                            |                 | (for long PCR products) |
| Human genomic DNA          | 5 - 500 ng      | (100 - 500 ng)          |
| <i>E. coli</i> genomic DNA | 100 pg - 200 ng | (10 - 200 ng)           |
| Plasmid DNA                | 10 pg - 10 ng   | (1 - 10 ng)             |
| cDNA                       | 25 - 750 ng     | (250 - 750 ng)          |

Note: Do not use templates containing uracil, such as bisulfite-treated DNA.

## Electrophoresis of amplified products:

TAE Buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using TaKaRa Ex Premier DNA Polymerase Dye plus. The vivid green dye separates into blue and yellow dye fronts when the PCR product is run on an agarose gel.

Note: TBE Buffer may result in DNA band patterns that are enlarged at the bottom of the gel.

## Termini of amplified products:

Most PCR products amplified with TaKaRa Ex Premier DNA Polymerase Dye plus have blunt-end termini.

Use the Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) for cloning into a blunt-end vector.

In-Fusion® cloning with In-Fusion Snap Assembly Master Mix series (Cat. #638947/638948/638949/638943/638944/638954/638955/638956) also can be recommended.

If you want to clone amplified products into T-vectors, dA addition to the 3' end is needed.

## Restriction enzyme digestion:

Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of TaKaRa Ex Premier DNA Polymerase Dye plus from the reaction mixture by phenol/chloroform extraction or by using the NucleoSpin Gel and PCR Clean-Up (Cat. #740609.10/.50/.250) kit. A 3'-overhang produced with restriction enzyme, such as *Pst* I, may be removed by the 3'  $\rightarrow$  5' exonuclease activity of TaKaRa Ex Premier DNA Polymerase Dye plus if residual polymerase remains in the restriction enzyme reaction.

## cDNA template:

The cDNA products synthesized by some reverse transcription reagents cannot be recommended as the template since the PCR solution may aggregate.

Please use the PrimeScript™ II 1st strand cDNA Synthesis Kit (Cat. #6210A/B) and PrimeScript IV 1st strand cDNA Synthesis Mix (Cat. #6215A/B) for cDNA synthesis.

In-Fusion is a registered trademark of Takara Bio USA, Inc.  
Ex Premier and PrimeScript are trademarks of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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