

Confident detection of low-frequency mutations in cell-free DNA using SMARTer ThruPLEX technology with unique molecular tags



Suvarna Gandlur^{1*}, Anthony Popkie¹, Emmanuel Kamberov¹, Matthew Carroll¹, Dmitry Goryunov, Konstantinos Charizanis³, Fang Sun¹, Maria Dinkelmann⁴, Kamran Shazand⁵, Shawn Quinn², John Langmore¹

¹Takara Bio USA, Inc., Mountain View, CA 94043, USA • ²Curio Genomics, Inc., Dexter MI 48130, USA • ³Gemphire Therapeutics, Livonia, MI 48152 • ⁴Zomedica, Ann Arbor, MI 48108 • ⁵PacGenomics, Agoura Hills, CA 91301
*Corresponding author: suvarna_gandlur@takarabio.com

Abstract

Liquid biopsies provide a noninvasive method to acquire the genetic information provided in cell-free DNA (cfDNA). Access to this genetic information through next-generation sequencing (NGS) identifies mutations and alterations that may play a role in cancer and other diseases. The key to identifying rare mutations is improved sequencing accuracy and the ability to distinguish between biological and PCR duplicates. SMARTer® ThruPLEX® Tag-seq was developed with the addition of unique molecular tags (UMTs) to improve sequencing accuracy by accounting for polymerase and sequencing errors and to increase confidence in rare allele identification. Libraries were prepared with SMARTer ThruPLEX Tag-seq using 10–30 ng of Horizon Discovery's Multiplex I cfDNA Reference Standard Set containing six single nucleotide variants (SNV) for four different genes (*EGFR*, *KRAS*, *NRAS*, *PIK3CA*) present at 0.5–5% allele frequency. The libraries were enriched with either a 110-kb or 240-kb custom panel or the Agilent ClearSeq Comprehensive Cancer Panel. Enriched libraries were sequenced with an average total read coverage of approximately 5,000X and analyzed with and without the UMTs. Deduplication without molecular tags reduced coverage to 295X; whereas, deduplication with UMTs allowed for a separation of biological duplicates from PCR duplicates and increased coverage to 2,110X, a significant reduction in false positives, 73% elimination of background noise, and a 10-fold increase in unique coverage compared to deduplication without UMTs. Employing UMT consensus reads, the sensitivity to detect 70 SNVs at 1% minor allele frequency (MAF) was increased from 30% to 95% reads and at 0.2% MAF, increased from 7% to 75% and false positive calls reduced by 32X. Therefore, use of UMTs in the preparation of NGS libraries from cfDNA enhances sequencing accuracy: by distinguishing between biological duplicates and PCR duplicates, increasing read coverage and decreasing background noise, reducing false positives, and in more confident mutation calls.

1 SMARTer ThruPLEX chemistry

SMARTer ThruPLEX Tag-seq library preparation is a simple, three-step, two-hour, one-tube reaction process with no intermediate cleanups. The repair and ligation reactions use proprietary stem-loop adaptors with added molecular tags that ligate only to the 5' ends of the cfDNA, leaving the 3' ends to be extended and indexed by high-fidelity amplification. (Figure 1).




Figure 1. SMARTer ThruPLEX chemistry. SMARTer ThruPLEX Tag-seq starts with cell-free DNA or fragmented double-stranded DNA. In the first step, DNA fragments are repaired in a highly efficient process. Following repair, stem-loop adapters are ligated to DNA fragments to attach the Illumina® adapter sequence and a unique molecular tag (UMT) to each side of the DNA fragment. In the final reaction, indexing primers containing Illumina P5 and P7 sequences are used to complete the library structure and amplify the library fragments.

2 Unique molecular tags provide confident variant detection

To establish the variant detection performance of the SMARTer ThruPLEX Tag-seq Kit, the limits of variant detection were measured using cfDNA reference standards engineered with variants at different allele frequencies. This effort involved the steps of library preparation, target enrichment, sequencing, and data analysis (Figure 2).



Figure 2. Complete workflow from sample to data analysis.

3 Uniform distribution of unique molecular tags uniquely label fragments




Figure 3. Incorporation of unique molecular tags. A. SMARTer ThruPLEX Tag-seq provides >16 million unique combinations of molecular tags. More than 4,000 UMTs are available on each side of the fragment. B. After sequencing and counting of the reads, more than 99% of the UMT counts are distributed within twofold of the mean, demonstrating unbiased incorporation.

4 Variant detection using Horizon standards

Six Horizon cfDNA reference standards were tested, ranging from 0% (wild type) to 5% MAF. All six variants were called at their expected frequencies with sensitivity and specificity over 99% (Figure 4). By combining deep sequencing with the SMARTer ThruPLEX Tag-seq Kit, it was possible to detect mutations present at 0.5% allele frequency using a starting input of just 10 ng of DNA. Lower detection limits can be achieved, depending on sample quality, input amount, capture efficiency, sequencing depth, and data processing algorithms.

DNA Input	Enrichment	Sequencing	Sample MAF	Detected variant frequency						Sensitivity	Specificity
				<i>EGFR L858R</i>	<i>EGFR T790M</i>	<i>KRAS G12D</i>	<i>NRAS A60T</i>	<i>NRAS Q61K</i>	<i>PIK3CA E545K</i>		
30 ng	110 kb panel, Agilent SureSelect	~1,000X mean unique coverage, NextSeq 500	5%	3.7%	5.7%	6.1%	4.4%	5.9%	5.6%	100.0%	99.8%
			1%	0.5%	1.4%	1.5%	1.7%	1.2%	1.0%	100.0%	99.9%
			WT**	0%*	0%*	0%*	0%*	0%*	0%*		
10 ng	240 kb panel, Roche NimbleGen	~500X mean unique coverage, HiSeq 2500	2.5%	2.3%	1.0%	2.5%	3.6%	2.3%	1.6%	100.0%	99.6%
			1%	1.4%	0.6%	1.3%	0.9%	0.4%	1.7%	100.0%	99.8%
			0.5%	1.4%	0.2%	0.9%	0.8%	1.3%	1.1%	100.0%	99.8%

*Not detected
**100% wild-type negative control

Figure 4. UMTs provide excellent variant detection. Horizon Multiplex I cfDNA Reference Standards (Horizon HD780) were used as is or titrated using the wild-type reference standard to generate samples at additional allele frequencies. Variants were detected at their expected mean allele frequencies (MAF) with high sensitivity and specificity.

5 Reduction of background noise




Figure 5. UMTs reduce background errors for more confident identification of variants. Panel A. Aligned raw reads show that the expected EGFR T790M mutation (gold dot) is obscured by false positive noise (gray dots), making it difficult to distinguish the true mutation from false positives. Panel B. In contrast, by using UMTs in the same sample, the consensus reads show a dramatic reduction in the level of background errors, and as a result, clearer separation of the true mutation from the noise.

6 UMTs improve signal-to-noise ratio




Figure 6. UMTs greatly increase signal-to-noise ratio. For each sample and processing method, the signal was calculated by taking the average of the allele frequency detected at six expected mutation positions, and the noise was calculated by averaging the allele frequency detected outside the six positions across the entire captured region. The results show a 3X to 6X improvement in signal-to-noise ratio when UMTs were utilized for error-correction during data processing.

7 Variant detection using cell-free DNA spike-in samples

Baylor Miraca Genetics laboratory prepared human cfDNA spike-in samples by combining plasma-derived cfDNA from two distinct donors. Libraries were prepared with SMARTer ThruPLEX Tag-seq using 30 ng of cfDNA isolated from plasma, enriched with a 250-kb custom panel and sequenced to 30,000X coverage. The data was analyzed using the Curio Genomics software, and 70 unique variants were identified.




Figure 7. UMTs increase unique coverage and sensitivity while reducing false positives. Panel A. SMARTer ThruPLEX Tag-seq libraries with UMTs show a 10-fold increase in unique coverage by utilizing biological duplicates (orange) compared to libraries analyzed without UMTs. Panel B. SMARTer ThruPLEX Tag-seq libraries with UMTs show higher sensitivity when detecting low-frequency alleles (<1%) compared to libraries without UMTs. Panel C. SMARTer ThruPLEX Tag-seq libraries with UMTs show higher specificity when compared to libraries without UMTs.

Conclusions

Equipped with more than 16 million UMTs, the SMARTer ThruPLEX Tag-seq library preparation kit is a powerful tool for confident detection of low-frequency variants. SMARTer ThruPLEX Tag-seq technology's highly efficient chemistry and single-tube workflow work together to preserve molecular complexity, allowing researchers to discover more from precious samples using just 1 to 50 ng of DNA. The kit provides the freedom to use any commercially available capture panels or to design custom capture panels to interrogate genomic regions of interest that span hundreds of genes. These panels also allow structural studies of the genomic regions. Lower detection limits can be achieved, depending on sample quality and input amount, capture efficiency, and sequencing depth.

Takara Bio USA, Inc.
United States/Canada: +1 800.662.2566 • Asia Pacific: +1 800.919.7300 • Europe: +33 (0)1.3604.8880 • Japan: +81 (0)77.565.6999
FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. © 2017 Takara Bio Inc. All Rights Reserved. All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com



To download a copy of this poster, visit:
<https://go.takarabio.com/ABRF-2018.html>

800.662.2566
Visit us at takarabio.com

Clontech Takara cellartis