

Advantage® Genomic LA Polymerase Mix Protocol-at-a-Glance

(PT3885-2)

Advantage Genomic LA Polymerase Mix (Cat. No. 639152) is composed of a full-length *Taq* DNA polymerase, a small amount of proofreading enzyme, and a hot-start antibody. This enzyme is suitable for hot-start PCR due to the presence of a neutralizing monoclonal antibody that binds and inhibits the DNA polymerase at ambient temperatures. Nonspecific amplification due to mispriming and/or formation of primer-dimers before thermal cycling is prevented because the antibody inhibits the polymerase activity by binding to the enzyme until the temperature rises.

The key features of this enzyme are synthesis of products up to 30 kb using human genomic DNA template and up to 48 kb using noncomplex templates; 6.5X higher fidelity than wild-type *Taq* DNA polymerase; and less optimization and greater product yields than other "long and accurate polymerases" due to the robust enzyme/buffer system.

The protocol for using the Advantage Genomic LA Polymerase Mix is outlined below.

PCR Protocol

1. Thaw all reagents on ice.
2. Assemble the PCR reactions by adding each of the components listed in Table I. This can be done at room temperature.

TABLE I. COMPONENTS OF PCR REACTIONS	
Reagents	Volume/Final Concentration per Reaction
Sterile deionized H ₂ O	up to 25 µl
10X Advantage Genomic LA Buffer (Mg ²⁺) ¹	2.5 µl
dNTP Mixture (10 mM each)	1 µl
Primer 1	0.2–1.0 µM
Primer 2	0.2–1.0 µM
Advantage Genomic LA Polymerase Mix (5 units/µl) ²	0.25 µl
Template DNA	0.1–1 µg

¹ 500 µl of Advantage Genomic LA Buffer is supplied in a vial. This buffer contains 25 mM MgCl₂.

² The enzyme is supplied in a storage buffer containing 20 mM Tris-HCl (pH 8.0), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Tween 20, 0.5% Nonidet P-40, and 50% glycerol.

3. Mix well and spin down briefly to collect all the liquid at the bottom of the wells.
4. Place the plates in a thermal cycler and begin thermal cycling immediately, using the following parameters (for amplification of a 17.5 kb fragment).
 - 94°C for 1 min
 - 30 cycles
 - 98°C for 10 sec
 - 68°C for 15 min
 - 72°C for 10 min

Notes:

- The denaturation conditions vary depending upon the thermal cycler and the type of plate used. We recommend 10–30 sec at 94°C or 1–10 sec at 98°C.
- The extension time should be determined by the size of the amplicon. As a guideline, use 1 min per 1 kb of amplicon.



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

(PR143852; published April 2011)

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Your use of this product is subject to compliance with any applicable licensing requirements described on the product's web page at www.clontech.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech, the Clontech logo, Advantage, and TaqStart are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2011 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.