

I. List of Components

- Terra Direct PCR Polymerase Mix (1.25 U/μl)
- 2X Terra Direct PCR Buffer (with Mg²⁺, dNTPs)—For Genotyping
- Terra Direct PCR—Tissue Extraction Buffer
- **NOTE:** Store at room temperature after thawing.
- Proteinase K.
- 5X Loading Dye—For Genotyping

II. Additional Materials Required

- PCR-grade water (Avoid using autoclaved H₂O; the recycled steam in some autoclaves can introduce contaminants that may interfere with PCR.)
- Gene-specific PCR primers (with T_m ≥ 60°C; avoid using primers that contain inosine)
- PCR reaction tubes or plate
- Aerosol-resistant pipette tips preferably equipped with hydrophobic filters.

III. General Considerations

Follow the recommendations in Option A if you are performing direct PCR with the whole-tissue samples and Option B if you are performing PCR with crude tissue extracts:

Option A

Direct PCR with Whole-Tissue Samples

Use the following as a guide to help you determine the appropriate amount of whole-tissue sample to *add directly to* your PCR reaction.

When using:

- mouse tail biopsies, add ≤ 1 mm of tissue.
- mouse ear biopsies, add ≤ 2 mm diameter disc of tissue.
- mouse organs, add ≤ 1.5 mm³.

Option B

Direct PCR with Crude Tissue Extracts

For tissue extracts, you may use any sample preparation method that is appropriate for your sample type.

We have found the following method works well for mouse tissue extracts:

- i. Add 20 μl of Terra Direct PCR—Tissue Extraction Buffer* to the mouse tissue.
- ii. Add 0.5 μl Proteinase K and mix by vortexing briefly.
- iii. Incubate the reaction for 5 min at room temp, and then for 2 min at 98°C.
- iv. Use ≤ 5 μl of the crude extract in the PCR reaction (see Table 1, next page).

* Perform the above procedure at room temperature. If a precipitate has formed in the Tissue Extraction Buffer, heat the tube at 60°C and invert gently until the precipitate has completely dissolved.

IV. Protocol

A. PCR Reaction Set-Up

1. Prepare the PCR reaction on ice by adding each component indicated in Table 1.
2. Briefly spin the tube in a microcentrifuge, and begin thermal cycling using the guidelines provided below.

Table 1. PCR Reaction Components

Reagent	Amount	Final Conc.
Sterile H ₂ O	Up to 50 µl	—
2X Terra Direct PCR Buffer—For Genotyping ^a	25 µl	1X ^b
Primer 1	15 pmol	0.3 µM
Primer 2	15 pmol	0.3 µM
Tissue Sample/Extract	≤ 5 µl ^c	
Terra Direct PCR Polymerase Mix	1 µl	1.25 U
Total volume	50 µl^d	

^a The buffer may turn white when left at room temperature; this is not a problem as it will become clear when placed on ice.

^b In the PCR reaction, the final concentration of Mg²⁺ is 2 mM and the final concentration of each dNTP is 400 mM.

^c See the ‘General Considerations’ section, above, for suggested amounts of different sample types.

^d For 25 µl reactions, be sure to add only half the amount indicated for each reagent.

B. Recommended Cycling Conditions

Use the following cycling conditions when setting up your initial experiments. These are general guidelines—the optimal conditions may vary.

3-Step PCR (recommended)

- 98°C for 2 min*
- 30–40 cycles:
 - 98°C for 10 sec
 - 60°C for 15 sec
 - 68°C for 1 min/kb

2-Step PCR

- 98°C for 2 min*
- 30–40 cycles:
 - 98°C for 10 sec
 - 68°C for 1 min/kb

* The initial denaturation step must be performed at 98°C for 2 min in order to denature the hot start antibody.

C. Post-PCR Considerations

1. When running your samples on agarose gels, **we recommend using the 5X Loading Dye—For Genotyping supplied in the kit.** This loading dye is formulated to provide improved sample resolution on agarose gels.
2. Use TAE running buffer when visualizing your PCR products by agarose gel electrophoresis. The use of TBE is not recommended as it causes spreading of the DNA bands toward the bottom of the gel.
3. PCR products produced by Terra Direct PCR contain 3' A-overhangs, making them compatible with T/A cloning.

Contact Us	
Customer Service/Ordering	Technical Support
tel: 800.662.2566 (toll-free)	tel: 800.662.2566 (toll-free)
fax: 800.424.1350 (toll-free)	fax: 800.424.1350 (toll-free)
web: takarabio.com/service	web: takarabio.com/support
e-mail: ordersUS@takarabio.com	e-mail: technical_support@takarabio.com

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.