

I. Introduction

Lyophilized PrimePath™ Probe qPCR Mix (Cat. No. 638351) enables the user to perform accurate probe-based qPCR on purified DNA samples using a lyophilized master mix. The cakes allow for flexibility in handling larger sample volumes, in addition to the convenience of room-temperature storage and transport.

II. Required Materials

This protocol applies to the following Takara Bio products:

- Lyophilized PrimePath Probe qPCR Mix* 3 x 8-tube strips

*Store at room temperature. Store opened, unused cakes in a desiccator.

Additional materials required:

- Primers and probes
- Micropipette and tips (with hydrophobic filters)
- Vortex mixer
- Benchtop centrifuge for tubes or plates
- 1.5 ml Eppendorf tubes, 200 µl PCR tubes, or 200 µl PCR plates for sample preparation
- Tubes or plates for real-time PCR with optical seals or caps
- A real-time PCR machine (e.g., CFX96 Real-Time PCR Detection System [Bio-Rad], QuantStudio 3 or 5 [Thermo Fisher Scientific])

III. Protocol

A. qPCR Reaction Mix Preparation

1. Prepare the reaction mix for the qPCR reaction as illustrated in the example below (optimization is recommended).

Example:

- **qPCR reaction mixture (1 rxn)**

1 cake	Lyophilized PrimePath Probe qPCR Mix
1 µl	20X primer/probe mix (provided by user)
17 µl	DNase/RNase-free H ₂ O
2 µl	DNA sample*
20 µl	Total volume

*The DNA range will be dependent on your individual assay and testing conditions. You will need to optimize the volume for your experimental parameters.

2. Vortex until the cake is completely dissolved then spin down.
3. Wait 5 min, then repeat Step 2.

NOT A SAFE STOPPING POINT: Proceed to the qPCR reaction protocol (Section B).

B. qPCR Reaction Protocol

1. Run the standard cycling protocol according to the instruction manual of your real-time qPCR instrument.

The example PCR program below was used for detecting a multiplex assay of five targets (*Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, a spike-in control, and Human RNaseP target) and high-quality results were obtained.

NOTE: Further optimization may be required for your experiment.

Example:

95°C	2 min	Initial denaturation
40 cycles:		
95°C	10 sec	Denaturation
60°C	30 sec	
Anneal, extend, and capture		

2. After the reaction is complete, check the amplification curve. Confirm that the analytical parameters are appropriate and that the Cq value has been calculated.

NOTE: If the default setting analysis does not work, perform manual analysis per the instrument’s instruction manual.

Contact Us	
Customer Service/Ordering	Technical Support
tel: 800.662.2566 (toll-free)	tel: 800.662.2566 (toll-free)
fax: 800.424.1350 (toll-free)	fax: 800.424.1350 (toll-free)
web: takarabio.com/service	web: takarabio.com/support
e-mail: ordersUS@takarabio.com	e-mail: technical_support@takarabio.com

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product’s web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.