



	End of DD-tag	AgeI	EcoRI	XhoI	SpeI	XbaI	NotI	BamHI
3121	AAA CCG GAA	ACC GGT	GAA TTC	CTC GAG	ACT AGT	TCT AGA	GCG GCC GCG	GAT CCC

pLVX-PTuner Green Vector Map and Multiple Cloning Site (MCS).

Description

pLVX-PTuner Green is a bicistronic, HIV-1-based, lentiviral expression vector that allows you to precisely regulate the amount of your protein of interest in virtually any mammalian cell type, including primary cells. The vector encodes a 12 kDa, FKBP (L106P) destabilization domain (DD; 1) that is expressed as an N-terminal tag on your protein of interest, causing rapid degradation of the fusion protein. Degradation of the DD-tagged protein can be prevented by the addition of Shield1 stabilizing ligand to the medium. Shield1 is a membrane permeable molecule that binds to the DD tag, 'shielding' the fusion protein from proteasomal degradation.

pLVX-PTuner Green allows the simultaneous expression of your DD-tagged protein of interest and the fluorescent protein ZsGreen1 from the same bicistronic mRNA transcript. ZsGreen1 is a human codon-optimized variant of the reef coral *Zoanthus sp.* green fluorescent protein, ZsGreen (2). Expression of the bicistronic transcript is driven by the constitutively active human cytomegalovirus immediate early promoter ($P_{CMV IE}$) located just upstream of the DD tag coding sequence. An encephalomyocarditis virus (EMCV) internal ribosome entry site (IRES), positioned between the MCS and the ZsGreen1 gene (see above), facilitates cap-independent translation of ZsGreen1 from an internal start site at the IRES/ZsGreen1 junction (3). Because ZsGreen1 is unaffected by the DD tag, it can be used as an indicator of transfection and transduction efficiency, as well as a marker for selection and cell sorting.

pLVX-PTuner Green contains all of the viral processing elements necessary for the production of replication-incompetent lentivirus, as well as elements to improve viral titer, transgene expression, and overall vector function. The woodchuck hepatitis virus posttranscriptional

(PR113813; published 17 January 2011)



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

regulatory element (WPRES) promotes RNA processing events and enhances nuclear export of viral and transgene RNA (4), leading to increased viral titers from packaging cells and enhanced expression of your gene of interest in target cells. In addition, the vector includes a Rev-response element (RRE), which further increases viral titers by enhancing the transport of unspliced viral RNA out of the nucleus (5). Finally, pLVX-PTuner Green also contains a central polypurine tract/central termination sequence element (cPPT/CTS). During target cell infection, this element creates a central DNA flap that increases nuclear import of the viral genome, resulting in improved vector integration and more efficient transduction (6). The vector also contains a pUC origin of replication and an *E. coli* ampicillin resistance gene (Amp^r) for propagation and selection in bacteria.

Use

pLVX-PTuner Green is available in the Lenti-X™ ProteoTuner™ Shield System N (w/ ZsGreen1) [Cat. No. 632175]. The vector is designed to constitutively coexpress a DD-tagged protein of interest and ZsGreen1 from $P_{CMV IE}$ when transduced into mammalian cells. In order to create your N-terminally tagged protein of interest, your gene of interest must be cloned into the MCS in the same reading frame as the DD tag sequence, and it must contain a stop codon at the end of its coding sequence. Before it can be transduced into target cells, the vector must be cotransfected into 293T cells with our Lenti-X™ HTX Packaging System (Cat. Nos. 631247 and 631249) and packaged into viral particles. This packaging system allows the safe production of high titer, infectious, replication-incompetent, VSV-G pseudotyped lentiviral particles that can infect a wide range of cell types, including non-dividing and primary cells (7).

ZsGreen1 is the brightest commercially available green fluorescent protein. The presence of this protein allows the selection of stable transductants by flow cytometry (or other detection methods) with standard FITC filter sets (ZsGreen1 has an excitation maximum of 493 nm and an emission maximum of 505 nm).

When cells expressing a DD-tagged protein of interest are grown in medium containing Shield1, the ligand binds to the DD tag and protects the fusion protein from degradation. As a result, the protein quickly accumulates inside the cells in amounts that are directly proportional to the concentration of Shield1 in the medium. If the cells are subsequently grown in medium lacking Shield1, the DD tag is no longer stabilized, and the fusion protein is rapidly degraded. Because the effects of Shield1 are concentration-dependent and reversible, it is possible to fine-tune the amount of fusion protein present in the cells simply by adjusting the concentration of the ligand in the medium (1).

Location of Features

- 5' LTR: 1–635
- PBS (primer binding site): 636–653
- Ψ (packaging signal): 685–822
- RRE (Rev-response element): 1303–1536
- cPPT/CTS (central polypurine tract/central termination sequence): 2028–2151
- $P_{CMV IE}$ (human cytomegalovirus immediate early promoter): 2185–2788
- DD (FKBP-L106P destabilization domain): 2806–3129
- MCS (multiple cloning site): 3131–3169
- IRES (encephalomyocarditis virus internal ribosome entry site): 3175–3749
- ZsGreen1 (*Zoanthus sp.* green fluorescent protein): 3750–4445
- WPRES (woodchuck posttranscriptional regulatory element): 4459–5050
- 3' LTR: 5254–5890
- pUC origin of replication: 6360–7030 (complementary)
- Amp^r (ampicillin resistance gene; β -lactamase): 7175–8171 (complementary)

Selection of Stable Transfectants

- Selectable marker: ZsGreen1.

Propagation in *E. coli*

- Suitable host strains: Stellar™ Competent Cells.
- Selectable marker: plasmid confers resistance to ampicillin (100 μ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: high

Notes:

The vector sequence was compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced. The viral supernatants produced by this lentiviral vector could contain potentially hazardous recombinant virus. Due caution must be exercised in the production and handling of recombinant lentivirus. Appropriate NIH, regional, and institutional guidelines apply.

References

1. Banaszynski, L. *et al.* (2006) *Cell* **126**(5):995–1004.
2. BD Living Colors pIRES2-ZsGreen1 Vector (January 2005) *Clontechniques* **XX**(1):8.
3. Jang, S. K. *et al.* (1988) *J. Virol.* **62**(8):2636–2643.
4. Zufferey, R. *et al.* (1999) *J. Virol.* **73**(4):2886–2892.
5. Cochrane, A. W. *et al.* (1990) *Proc. Natl. Acad. Sci. USA* **87**(3):1198–1202.
6. Zennou, V. *et al.* (2000) *Cell* **101**(2):173–185.
7. Wu, X. *et al.* (2000) *Mol. Ther.* **2**(1):47–55.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

The RCFPs (including DsRed-Express and DsRed-Express2) are covered by one or more of the following U.S. Patent Nos. 7,166,444; 7,157,565; 7,217,789; 7,338,784; 7,338,783; 7,537,915 6,969,597, 7,150,979 and 7,442,522.

Clontech has a license to sell products containing WPRE, under the terms described below. Any use of WPRE outside of Clontech's product or the product's intended use, requires a license as detailed below. Before using the product containing WPRE, please read the following license agreement. If you do not agree to be bound by its terms, contact Clontech within 10 days for authorization to return the unused product containing WPRE and to receive a full credit. Patents: The WPRE technology is covered by one or more of the following U.S. Patents and corresponding patent claims outside the U.S.: 6,136,597; 6,284,469; 6,312,912; 6,287,814, issued to The Salk Institute for Biological Studies; Individual License Agreement: Clontech grants you a non-exclusive license to use the enclosed product containing WPRE in its entirety for its intended use. The product is being transferred to you in furtherance of, and reliance on, such license. Any use of WPRE outside of Clontech's product or the product's intended use, requires a license from the Salk Institute for Biological Studies. Termination of License: This license agreement is effective until terminated. You may terminate it at any time by destroying all products containing WPRE in your control. It will also terminate automatically if you fail to comply with the terms and conditions of the license agreement. You shall, upon termination of the license agreement, destroy all products containing WPRE in your control, and so notify Clontech in writing. This License shall be governed in its interpretation and enforcement by the laws of the State of California. Contact for WPRE Licensing: The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: 858.453.4100 ext. 1275, Fax: 858.546.8093

This product and its use are the subject of U.S. Pat. No. 6,682,907. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot disclose information, sell or otherwise transfer this product, its components or materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for any commercial purposes. If the buyer is not willing to accept the limitations of this limited use statement, Clontech is willing to accept return of the product with a full refund. For information on purchasing a license to the DNA-Flap technology for purposes other than research, contact the Transfer of Technology Office, Institut Pasteur, 28 rue du Docteur Roux, 75 724 Paris Cedex 15 (www.pasteur.fr).

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com.

For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2011 Clontech Laboratories, Inc.