



Additionally, the viral genomic transcript contains the necessary viral RNA processing elements, including the LTRs, packaging signal ( $\Psi^+$ ), and tRNA primer-binding site. pRetroQ-mCherry-C1 contains a bacterial origin of replication, an *E. coli* Amp<sup>r</sup> gene for propagation and selection in bacteria, and an SV40 origin for replication in mammalian cells expressing the SV40 large T antigen.

### Use

pRetroQ-mCherry-C1 is designed to efficiently deliver and express C-terminal mCherry fusions into primary cells or cells that are difficult to transfect. C-terminal mCherry fusion proteins retain the fluorescent properties of the native protein, allowing the *in vivo* localization of the fusion protein. The gene of interest should be cloned into pRetroQ-mCherry-C1 so that it is in frame with the mCherry coding sequence. The inserted sequence does not require an initiation codon (ATG) or a stop codon (TAA, TAG, TGA); however, if you don't want to use the stop codons downstream of the MCS (see map), you can add a stop codon to the end of your gene of interest. The recombinant mCherry vector can be transduced or transfected into mammalian cells. If required, stable transformants can be selected using puromycin. pRetroQ-mCherry-C1 can also be used simply to express mCherry in a cell line of interest (e.g., as an infection marker).

Before pRetroQ-mCherry-C1 can be transduced into mammalian cells, it must be transfected into a packaging cell line (such as the RetroPack™ PT67 Cell line (Cat. No. 631510), AmphiPack™ -293 (Cat. No. 631505), EcoPack™ 2-293 (Cat. No. 631507), Pantropic Expression System (Cat. No. 631512), or Retro-X™ Universal Packaging System (Cat. No. 631530)). The packaging cell line supplies the viral structural genes (gag, pol, and env) necessary for particle formation and replication that pRetroQ-mCherry-C1 lacks, allowing RNA from the vector to be packaged into non-infectious, replication-incompetent retroviral particles. Once a high-titer supernatant is produced, these retroviral particles can infect target cells and transmit the gene of interest, but they cannot replicate within the target cells due to the absence of viral structural genes. The separate introduction and integration of the structural genes into the packaging cell line minimizes the chances of producing replication-competent virus caused by recombination events during cell proliferation.

### Location of features

- 5' LTR (CMV/MSV): 1–728
- $\Psi^+$  (extended packaging signal): 757–1566
- $P_{CMVIE}$  (human cytomegalovirus immediate-early promoter ): 1582–2170
- mCherry fluorescent protein gene: 2194–2901
  - Start codon: (ATG): 2194–2196
  - Last codon: 2899–2901
- MCS (multiple cloning site): 2911–2967
- $P_{PGK}$  (PGK promoter): 2991–3499
- Puro<sup>r</sup> (puromycin resistance gene ): 3520–4119
- 3' LTR (MMLV; deletion in U3): 4304–4736
  - PolyA signal: 4562–4577
- $P_{SV40}$  (SV40 promoter): 5016–5283
- SV40 origin of replication: 5237–5302
- ColE1 origin of replication: 5623
- Amp<sup>r</sup> (ampicillin resistance gene;  $\beta$ -lactamase): 6383–7243 (complementary)

### Propagation in *E. coli*

- Suitable host strains: DH5 $\alpha$ ™, Fusion Blue, and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100  $\mu$ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: ColE1
- Copy number: low

### Excitation and emission maxima of mCherry

- Excitation maximum = 587 nm
- Emission maximum = 610 nm

**Notes:**

The vector sequence was compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

The viral supernatants produced by this retroviral vector could, depending on your cloned insert, contain potentially hazardous recombinant virus. Due caution must be exercised in the production and handling of recombinant retrovirus. Appropriate NIH, regional, and institutional guidelines apply.

**References**

1. Shaner, N.C. *et al.* (2004) *Nat. Biotechnol.* **22**(12):1567–1572.
2. Gorman, C. (1985). In DNA cloning: A practical approach, vol. II. Ed. D.M. Glover. (IRL Press, Oxford, U.K.) pp. 143–190.
3. Kinsella, T. M. & Nolan G. P. (1996) *Hum. Gene Ther.* **7**(12):1405–1413.
4. Ory, D. S. *et al.* (1996) *Proc. Nat. Acad. Sci. USA* **93**(21):11400–11406.
5. Pear, W. S. *et al.* (1993) *Proc. Natl. Acad. Sci. USA* **90**(18):8392–8396.
6. Yang, S. *et al.* (1999) *Hum. Gene Ther.* **10**(1):123–132.
7. Emerman, M. & Temin, H. M. (1984) *Cell* **39**(3):449–467.

**Notice to Purchaser**

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

For Research Use Only, Not for Use in Humans.

DH5 $\alpha$ <sup>TM</sup> is a trademark of Invitrogen Corporation.

**DsRed-Monomer and Fruit Fluorescent Proteins:**

The DsRed-Monomer and the Fruit Fluorescent Proteins are covered by one or more of the following U.S. Patents: 7,005,511; 7,157,566; 7,393,923 and 7,250,298.

**Living Colors<sup>®</sup> Fluorescent Protein Products:**

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com) or [click here for more information](#).

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2009 Clontech Laboratories, Inc.