



# Safety Data Sheet

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision Date 2025-01-16

Revision Number 6

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

**Product Code** ST1843  
**Product Name** Cell Extraction Enzyme  
**Pure substance/mixture** Mixture  
**1.2. Relevant identified uses of the substance or mixture and uses advised against**

**Identified uses** For research use only. Not for use in diagnostic procedures  
**Uses advised against** No information available

### 1.3. Details of the supplier of the safety data sheet

#### Supplier

##### USA:

Takara Bio USA, Inc.  
2560 Orchard Parkway  
San Jose, CA 95131, USA  
Phone: 800.662.2566/888.251.6618  
Web: [www.takarabio.com](http://www.takarabio.com)

##### Europe:

Takara Bio Europe S.A.S.  
34, Rue de la Croix de Fer  
78100 Saint-Germain-en-Laye, France  
Phone: +33.1.39.04.68.80  
Web: [www.takarabio.com](http://www.takarabio.com)

##### Europe:

Takara Bio Europe AB  
Arvid Wallgrens Backe 20,  
SE-413 46 Göteborg, Sweden  
Phone: +46.31.758.09.00  
Web: [www.takarabio.com](http://www.takarabio.com)

##### India:

DSS Takara Bio India Pvt. Ltd.  
A-5 Mohan Co-operative Industrial Estate, Mathura Road,  
New Delhi 110044, India  
Phone: +91.1800.212.4922 (Toll free)  
Web: [www.takarabio.com](http://www.takarabio.com)

For further information, please contact:

### 1.4. Emergency telephone number

**Emergency telephone** In case of emergency, call PERS (Professional Emergency Resource Services)  
1-800-633-8253 (US) or 801-629-0667 (international).

Italy	Marco Marano
-------	--------------

	CAV "Osp. Pediatrico Bambino Gesù" Dip. Emergenza e Accettazione DEA Roma, Piazza Sant'Onofrio,4 00165 0668593726
--	---

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to

Regulation (EC) No. 1272/2008 [CLP]

This mixture is classified as not hazardous according to regulation (EC) 1272/2008 [CLP]

### 2.2. Label elements

This mixture is classified as not hazardous according to regulation (EC) 1272/2008 [CLP]

#### Hazard statements

This mixture is classified as not hazardous according to regulation (EC) 1272/2008 [CLP]

EUH210 - Safety data sheet available on request

### 2.3. Other hazards

No information available.

**Endocrine Disruptor Information** Contains a known or suspected endocrine disruptor.

Chemical name	EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances of Very High Concern (SVHC) for Authorisation	EU - REACH (1907/2006) - Endocrine Disruptor Assessment List of Substances
Poly(oxyethylene) octylphenyl ether (Triton X-100)	Endocrine disrupting properties	-

Chemical name	Endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100(3) or Commission Regulation (EU) 2018/605(4)
Poly(oxyethylene) octylphenyl ether (Triton X-100)	Endocrine disrupting properties

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

Not applicable

### 3.2 Mixtures

Chemical name	Weight-%	REACH registration number	EC No (EU Index No)	Classification according to Regulation (EC) No. 1272/2008 [CLP]	Specific concentration limit (SCL)	M-Factor	M-Factor (long-term)
Glycerol 56-81-5	20 - 30	No data available	200-289-5	No data available	-	-	-
Calcium chloride 10043-52-4	0.1 - 1	No data available	233-140-8 (017-013-00-2)	Eye Irrit. 2 (H319)	-	-	-
Poly(oxyethylene) octylphenyl ether (Triton X-100) 9002-93-1	0.1 - 1	No data available	-	No data available	-	-	-

**Full text of H- and EUH-phrases: see section 16**

#### Acute Toxicity Estimate

If LD50/LC50 data is not available or does not correspond to the classification category, then the appropriate conversion value from CLP Annex I, Table 3.1.2, is used to calculate the acute toxicity estimate (ATEmix) for classifying a mixture

based on its components

Chemical name	Oral LD50 mg/kg	Dermal LD50 mg/kg	Inhalation LC50 - 4 hour - dust/mist - mg/L	Inhalation LC50 - 4 hour - vapor - mg/L	Inhalation LC50 - 4 hour - gas - ppm
Glycerol 56-81-5	27200	10000	5.85	No data available	No data available
Calcium chloride 10043-52-4	1000	5000	No data available	No data available	No data available
Poly(oxyethylene) octylphenyl ether (Triton X-100) 9002-93-1	1800	No data available	No data available	No data available	No data available

This product contains one or more candidate substance(s) of very high concern (Regulation (EC) No. 1907/2006 (REACH), Article 59)

Chemical name	CAS No.	SVHC candidates
Poly(oxyethylene) octylphenyl ether (Triton X-100)	9002-93-1	X

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

<b>Inhalation</b>	Remove to fresh air.
<b>Eye contact</b>	Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.
<b>Skin contact</b>	Wash skin with soap and water. In the case of skin irritation or allergic reactions see a physician.
<b>Ingestion</b>	Rinse mouth.

### 4.2. Most important symptoms and effects, both acute and delayed

<b>Symptoms</b>	No information available.
-----------------	---------------------------

### 4.3. Indication of any immediate medical attention and special treatment needed

<b>Note to physicians</b>	Treat symptomatically.
---------------------------	------------------------

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

<b>Suitable Extinguishing Media</b>	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
<b>Large Fire</b>	CAUTION: Use of water spray when fighting fire may be inefficient.
<b>Unsuitable extinguishing media</b>	Do not scatter spilled material with high pressure water streams.

### 5.2. Special hazards arising from the substance or mixture

<b>Specific hazards arising from the</b>	No information available.
--	---------------------------

chemical

### 5.3. Advice for firefighters

**Special protective equipment and precautions for fire-fighters** Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

**Personal precautions** Ensure adequate ventilation.

**For emergency responders** Use personal protection recommended in Section 8.

### 6.2. Environmental precautions

**Environmental precautions** See Section 12 for additional Ecological Information.

### 6.3. Methods and material for containment and cleaning up

**Methods for containment** Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Take up mechanically, placing in appropriate containers for disposal.

**Prevention of secondary hazards** Clean contaminated objects and areas thoroughly observing environmental regulations.

### 6.4. Reference to other sections

**Reference to other sections** See section 8 for more information. See section 13 for more information.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

**Advice on safe handling** Ensure adequate ventilation.

**General hygiene considerations** Handle in accordance with good industrial hygiene and safety practice.

### 7.2. Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep container tightly closed in a dry and well-ventilated place.

### 7.3. Specific end use(s)

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Exposure Limits

Chemical name	European Union	Austria	Belgium	Bulgaria	Croatia
Glycerol 56-81-5	-	-	TWA: 10 mg/m <sup>3</sup>	-	TWA: 10 mg/m <sup>3</sup>
Chemical name	Cyprus	Czech Republic	Denmark	Estonia	Finland

Glycerol 56-81-5	-	TWA: 10 mg/m <sup>3</sup> Ceiling: 15 mg/m <sup>3</sup>	-	TWA: 10 mg/m <sup>3</sup>	TWA: 20 mg/m <sup>3</sup>
Calcium chloride 10043-52-4	-	TWA: 5 mg/m <sup>3</sup> Ceiling: 4 mg/m <sup>3</sup>	-	-	-
Chemical name	France	Germany TRGS	Germany DFG	Greece	Hungary
Glycerol 56-81-5	TWA: 10 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup> Peak: 400 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup>	-
Chemical name	Ireland	Italy MDLPS	Italy AIDII	Latvia	Lithuania
Calcium chloride 10043-52-4	-	-	-	TWA: 2 mg/m <sup>3</sup>	-
Chemical name	Luxembourg	Malta	Netherlands	Norway	Poland
Glycerol 56-81-5	-	-	-	-	TWA: 10 mg/m <sup>3</sup>
Chemical name	Portugal	Romania	Slovakia	Slovenia	Spain
Glycerol 56-81-5	TWA: 10 mg/m <sup>3</sup>	-	TWA: 10 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup> STEL: 400 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup>
Chemical name	Sweden		Switzerland		United Kingdom
Glycerol 56-81-5	-		TWA: 50 mg/m <sup>3</sup> STEL: 100 mg/m <sup>3</sup>		TWA: 10 mg/m <sup>3</sup> STEL: 30 mg/m <sup>3</sup>

**Biological occupational exposure limits**

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

**Derived No Effect Level (DNEL)** No information available.

**Predicted No Effect Concentration (PNEC)** No information available.

**8.2. Exposure controls****Personal Protective Equipment**

**Eye/face protection** No special protective equipment required.

**Skin and body protection** No special protective equipment required.

**Respiratory protection** No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

**General hygiene considerations** Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls** No information available.

## SECTION 9: Physical and chemical properties

**9.1. Information on basic physical and chemical properties**

**Physical state** Liquid  
**Appearance** Clear, colorless  
**Color** No information available  
**Odor** Odorless  
**Odor Threshold** No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
<b>Melting point / freezing point</b>	No data available	None known
<b>Boiling point/boiling range (°C)</b>	No data available	None known

Flammability (solid, gas)	No data available	None known
Flammability Limit in Air		None known
Upper flammability limit:	No data available	
Lower flammability limit:	No data available	
Flash point	No data available	Open cup
Autoignition temperature	No data available	None known
Decomposition temperature		None known
pH	No data available	None known
pH (as aqueous solution)	No data available	No information available
Kinematic viscosity	No data available	None known
Dynamic Viscosity	No data available	None known
Water solubility	No data available	None known
Solubility in other solvents	No data available	None known
Partition coefficient	No data available	None known
Vapor pressure	No data available	None known
Relative density	No data available	None known
Bulk Density	No data available	
Liquid Density	No data available	
Vapor density	No data available	None known
Particle characteristics		
Particle Size	No information available	
Particle Size Distribution	No information available	

## 9.2. Other information

9.2.1. Information with regard to physical hazard classes  
Not applicable

9.2.2. Other safety characteristics  
No information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

Reactivity No information available.

### 10.2. Chemical stability

Stability Stable under normal conditions.

#### Explosion Data

Sensitivity to mechanical impact None.

Sensitivity to static discharge None.

### 10.3. Possibility of hazardous reactions

Possibility of hazardous reactions None under normal processing.

### 10.4. Conditions to avoid

Conditions to avoid None known based on information supplied.

### 10.5. Incompatible materials

Incompatible materials None known based on information supplied.

### 10.6. Hazardous decomposition products

Hazardous Decomposition Products None known based on information supplied.

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Information on likely routes of exposure

##### Product Information

Inhalation	Specific test data for the substance or mixture is not available.
Eye contact	Specific test data for the substance or mixture is not available.
Skin contact	Specific test data for the substance or mixture is not available.
Ingestion	Specific test data for the substance or mixture is not available.

#### Symptoms related to the physical, chemical and toxicological characteristics

Symptoms	No information available.
----------	---------------------------

#### Numerical measures of toxicity

##### Acute toxicity

The following values are calculated based on chapter 3.1 of the GHS document

ATEmix (oral)	38,501.00 mg/kg
ATEmix (dermal)	32,074.40 mg/kg
ATEmix (inhalation-dust/mist)	10.90 mg/l

##### Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Glycerol	= 27200 mg/kg ( Rat )	> 10 g/kg ( Rabbit )	> 5.85 mg/L ( Rat ) 4 h
Calcium chloride	= 1000 mg/kg ( Rat )	> 5000 mg/kg ( Rabbit )	-
Poly(oxyethylene) octylphenyl ether (Triton X-100)	= 1800 mg/kg ( Rat )	-	-

#### Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation	No information available.
Serious eye damage/eye irritation	No information available.
Respiratory or skin sensitization	No information available.
Germ cell mutagenicity	No information available.
Carcinogenicity	No information available.
Reproductive toxicity	No information available.

STOT - single exposure No information available.

STOT - repeated exposure No information available.

Aspiration hazard No information available.

## 11.2. Information on other hazards

### 11.2.1. Endocrine disrupting properties

Endocrine disrupting properties

### 11.2.2. Other information

Other adverse effects No information available.

## SECTION 12: Ecological information

### 12.1. Toxicity

#### Ecotoxicity

Unknown aquatic toxicity Contains 0.34423 % of components with unknown hazards to the aquatic environment.

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Glycerol	-	LC50: 51 - 57mL/L (96h, <i>Oncorhynchus mykiss</i> )	-	-
Calcium chloride	-	LC50: =10650mg/L (96h, <i>Lepomis macrochirus</i> )	-	LC50: 2280000 - 3948000µg/L (48h, <i>Daphnia magna</i> )

### 12.2. Persistence and degradability

Persistence and degradability No information available.

### 12.3. Bioaccumulative potential

Bioaccumulation There is no data for this product.

#### Component Information

Chemical name	Partition coefficient
Glycerol	-1.75

### 12.4. Mobility in soil

Mobility in soil No information available.

### 12.5. Results of PBT and vPvB assessment

PBT and vPvB assessment No information available.

Chemical name	PBT and vPvB assessment
Glycerol	The substance is not PBT / vPvB
Calcium chloride	The substance is not PBT / vPvB

### 12.6. Endocrine disrupting properties



**Endocrine disrupting properties** No information available.

**12.7. Other adverse effects**

No information available.

## SECTION 13: Disposal considerations

**13.1. Waste treatment methods**

**Waste from residues/unused products** Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

**Contaminated packaging** Do not reuse empty containers.

## SECTION 14: Transport information

**IATA**

14.1 UN number or ID number	Not regulated
14.2 UN proper shipping name	No information available
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable
14.6 Special precautions for user	
Special Provisions	None

**IMDG**

14.1 UN number or ID number	Not regulated
14.2 UN proper shipping name	No information available
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable
14.6 Special precautions for user	
Special Provisions	None
14.7 Maritime transport in bulk according to IMO instruments	No information available

**RID**

14.1 UN number or ID number	Not regulated
14.2 UN proper shipping name	No information available
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable
14.6 Special precautions for user	
Special Provisions	None

**ADR**

14.1 UN number or ID number	Not regulated
14.2 UN proper shipping name	No information available
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable
14.6 Special precautions for user	
Special Provisions	None

## SECTION 15: Regulatory information

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****European Union**

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work.

**Authorizations and/or restrictions on use:**

This product contains one or more substance(s) subject to authorization (Regulation (EC) No. 1907/2006 (REACH), Annex XIV)

This product contains one or more substance(s) subject to restriction (Regulation (EC) No. 1907/2006 (REACH), Annex XVII)

Chemical name	Restricted substance per REACH Annex XVII	Substance subject to authorization per REACH Annex XIV
Calcium chloride - 10043-52-4	75	-
Poly(oxyethylene) octylphenyl ether (Triton X-100) - 9002-93-1	-	42

**Persistent Organic Pollutants**

Not applicable

**Ozone-depleting substances (ODS) Regulation (EU) 2024/590**

Not applicable

**International Inventories**

TSCA	-
DSL/NDL	-
EINECS/ELINCS	-
ENCS	-
IECSC	-
KECI	-
PICCS	-
AICS	-

**Legend:**

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**EINECS/ELINCS** - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

**IECSC** - China Inventory of Existing Chemical Substances

**KECI** - Korean Existing Chemicals Inventory

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**15.2. Chemical safety assessment**

**Chemical Safety Assessment** No information available

**SECTION 16: Other information****Key or legend to abbreviations and acronyms used in the safety data sheet**

**Legend**

SVHC: Substances of Very High Concern for Authorization:

**Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION**

TWA	Time weighted average	STEL	Short term exposure limit
Ceiling	Maximum limit value	*	Skin designation
**	Hazard Designation	+	Sensitizers

Classification procedure	
Classification according to Regulation (EC) No. 1272/2008 [CLP]	Method Used
Acute oral toxicity	Calculation method
Acute dermal toxicity	Calculation method
Acute inhalation toxicity - gas	Calculation method
Acute inhalation toxicity - vapor	Calculation method
Acute inhalation toxicity - dust/mist	Calculation method
Skin corrosion/irritation	Calculation method
Serious eye damage/eye irritation	Calculation method
Respiratory sensitization	Calculation method
Skin sensitization	Calculation method
Mutagenicity	Calculation method
Carcinogenicity	Calculation method
Reproductive toxicity	Calculation method
STOT - single exposure	Calculation method
STOT - repeated exposure	Calculation method
Acute aquatic toxicity	Calculation method
Chronic aquatic toxicity	Calculation method
Aspiration hazard	Calculation method
Ozone	Calculation method

**Key literature references and sources for data used to compile the SDS**

Agency for Toxic Substances and Disease Registry (ATSDR)  
 U.S. Environmental Protection Agency ChemView Database  
 European Food Safety Authority (EFSA)  
 Environmental Protection Agency  
 Acute Exposure Guideline Level(s) (AEGL(s))  
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act  
 U.S. Environmental Protection Agency High Production Volume Chemicals  
 Food Research Journal  
 Hazardous Substance Database  
 International Uniform Chemical Information Database (IUCLID)  
 National Institute of Technology and Evaluation (NITE)  
 Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)  
 NIOSH (National Institute for Occupational Safety and Health)  
 National Library of Medicine's ChemID Plus (NLM CIP)  
 National Library of Medicine's PubMed database (NLM PUBMED)  
 U.S. National Toxicology Program (NTP)  
 New Zealand's Chemical Classification and Information Database (CCID)  
 Organization for Economic Co-operation and Development Environment, Health, and Safety Publications  
 Organization for Economic Co-operation and Development High Production Volume Chemicals Program  
 Organization for Economic Co-operation and Development Screening Information Data Set  
 World Health Organization

Revision Date 2025-01-16

**This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006****Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless

specified in the text.

**End of Safety Data Sheet**