Safety Data Sheet



Revision Date 2023-12-28 Revision Number 7

1. Identification

Product identifier

Product Name SMART-Seq ICELL8 Forward Indexing Primer Set 1-72

Other means of identification

Product Code ST1849

Synonyms No information available

Recommended use of the chemical and restrictions on use

Identified uses For research use only. Not for use in diagnostic procedures

Restrictions on use No information available

Details of the supplier of the safety data sheet

<u>Supplier</u>

USA: Takara Bio USA, Inc. 2560 Orchard Parkway San Jose, CA 95131, USA

Phone: 800.662.2566/888.251.6618

Web: www.takarabio.com

Emergency telephone number

Emergency telephone In case of emergency, call PERS (Professional Emergency Resource Services)

1-800-633-8253 (US) or 801-629-0667 (international).

2. Hazard(s) identification

Product Classification Data

| Acute toxicity - Inhalation (Dusts/Mists) | Category 4 |
|---|-------------|
| Serious eye damage/eye irritation | Category 2 |
| Carcinogenicity | Category 1B |
| Reproductive toxicity | Category 1B |

Hazards not otherwise classified (HNOC)

Not applicable

Label elements

Danger

Hazard statements

Harmful if inhaled

Causes serious eye irritation

May cause cancer

May damage fertility or the unborn child



Appearance Clear, colorless Physical state Liquid Odor Amine

Precautionary Statements - Prevention

Obtain special instructions before use

Do not handle until all safety precautions have been read and understood

Wear protective gloves/clothing and eye/face protection

Avoid breathing dust/fume/gas/mist/vapors/spray

Use only outdoors or in a well-ventilated area

Wash face, hands and any exposed skin thoroughly after handling

Precautionary Statements - Response

IF exposed or concerned: Get medical advice/attention

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

If eye irritation persists: Get medical advice/attention

IF INHALED: Remove person to fresh air and keep comfortable for breathing

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Other information

May be harmful in contact with skin.

3. Composition/information on ingredients

Substance

Not applicable.

<u>Mixture</u>

| Chemical name | CAS No. | Weight-% | Trade secret |
|-----------------------|---------|----------|--------------|
| N,N-Dimethylformamide | 68-12-2 | 40 - 50 | * |

4. First-aid measures

Description of first aid measures

General advice IF exposed or concerned: Get medical advice/attention. Show this safety data sheet to the

doctor in attendance.

Inhalation Remove to fresh air. If symptoms persist, call a physician. If breathing has stopped, give

artificial respiration. Get medical attention immediately.

Eye contact Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. Keep eye wide open while rinsing. Do not rub affected area. Get medical attention if irritation develops and

persists.

Skin contact Wash skin with soap and water.

Ingestion Rinse mouth. Never give anything by mouth to an unconscious person. Do NOT induce

vomiting. Get medical attention.

Self-protection of the first aider

Avoid contact with skin, eyes or clothing. Ensure that medical personnel are aware of the

material(s) involved, take precautions to protect themselves and prevent spread of contamination. Avoid breathing vapors or mists. Use personal protective equipment as

required. See section 8 for more information.

Most important symptoms and effects, both acute and delayed

Symptoms May cause redness and tearing of the eyes. Burning sensation. Coughing and/ or wheezing.

Difficulty in breathing.

Indication of any immediate medical attention and special treatment needed

Note to physiciansTreat symptomatically.

5. Fire-fighting measures

surrounding environment.

Large Fire CAUTION: Use of water spray when fighting fire may be inefficient.

Unsuitable extinguishing media Do not scatter spilled material with high pressure water streams.

Specific hazards arising from the

chemical

No information available.

Explosion Data

Sensitivity to mechanical impact None.

Sensitivity to static discharge No

None.

Special protective equipment and precautions for fire-fighters

Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear.

Use personal protection equipment.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Personal precautions Avoid contact with skin, eyes or clothing. Use personal protective equipment as required.

Ensure adequate ventilation. Avoid breathing vapors or mists.

Other information Refer to protective measures listed in Sections 7 and 8.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Pick up and transfer to properly labeled containers.

7. Handling and storage

Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Avoid contact with

skin, eyes or clothing. Do not eat, drink or smoke when using this product. Remove contaminated clothing and shoes. Avoid breathing vapors or mists. Ensure adequate ventilation. In case of insufficient ventilation, wear suitable respiratory equipment.

Conditions for safe storage, including any incompatibilities

Storage Conditions Store locked up. Keep containers tightly closed in a dry, cool and well-ventilated place.

Keep out of the reach of children.

8. Exposure controls/personal protection

Control parameters

Exposure Limits

| Chemical name | ACGIH TLV | OSHA PEL | NIOSH |
|-----------------------|------------|-------------------------------------|---------------------------|
| N,N-Dimethylformamide | TWA: 5 ppm | TWA: 10 ppm | IDLH: 500 ppm |
| 68-12-2 | Sk* | TWA: 30 mg/m ³ | TWA: 10 ppm |
| | | (vacated) TWA: 10 ppm | TWA: 30 mg/m ³ |
| | | (vacated) TWA: 30 mg/m ³ | - |
| | | (vacated) Sk* | |
| | | Sk* | |

Biological occupational exposure limits

| Chemical name | ACGIH |
|-----------------------|--|
| N,N-Dimethylformamide | 30 mg/L - urine (total N-Methylformamide) - end of shift |
| 68-12-2 | 30 mg/L - urine (N-Acetyl-S-(N-methylcarbamoyl)cysteine) |
| | - end of shift at end of workweek |

Appropriate engineering controls

Engineering controls Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection If splashes are likely to occur, wear safety glasses with side-shields.

Hand protection Wear suitable gloves.

Skin and body protectionWear suitable protective clothing.

exceeded or irritation is experienced, ventilation and evacuation may be required.

General hygiene considerations

Do not eat, drink or smoke when using this product. Wash hands before breaks and immediately after handling the product. Avoid contact with skin, eyes or clothing. Wear suitable gloves and eye/face protection.

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state Liquid

Appearance Clear, colorless

Color No information available

Odor Amine

Odor Threshold No information available

<u>Property</u> <u>Values</u> <u>Remarks • Method</u>

No data available рΗ None known Melting point / freezing point No data available None known Boiling point/boiling range (°C) No data available None known Flash point No data available Open cup **Evaporation Rate** No data available None known Flammability (solid, gas) No data available None known Flammability Limit in Air None known

Upper flammability limit:No data availableLower flammability limit:No data available

Vapor pressure No data available None known Vapor density No data available None known Relative density No data available None known Water solubility No data available None known Solubility in other solvents No data available None known Partition coefficient No data available None known Autoignition temperature No data available None known **Decomposition temperature** None known

Kinematic viscosity

No data available

None known

No data available

None known

Other information

Explosive properties

Oxidizing properties

No information available
VOC content
No information available
Liquid Density
No information available
Bulk Density
No information available

10. Stability and reactivity

Reactivity No information available.

Chemical stability Stable under normal conditions.

Possibility of hazardous reactions None under normal processing.

Conditions to Avoid Excessive heat.

Incompatible materialsNone known based on information supplied.

Hazardous decomposition products None known based on information supplied.

11. Toxicological information

Information on likely routes of exposure

Product Information

Inhalation Specific test data for the substance or mixture is not available. May cause irritation of

respiratory tract. Harmful by inhalation. (based on components).

Eye contact Specific test data for the substance or mixture is not available. Causes serious eye irritation.

(based on components). May cause redness, itching, and pain.

Skin contact Specific test data for the substance or mixture is not available. May cause irritation.

Prolonged contact may cause redness and irritation. May be harmful in contact with skin.

Ingestion Specific test data for the substance or mixture is not available. Ingestion may cause

gastrointestinal irritation, nausea, vomiting and diarrhea.

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms May cause redness and tearing of the eyes. Coughing and/ or wheezing.

Acute toxicity

Numerical measures of toxicity

The following values are calculated based on chapter 3.1 of the GHS document

ATEmix (oral) 5,833.30 mg/kg ATEmix (dermal) 2,291.70 mg/kg ATEmix (inhalation-dust/mist) 3.12 mg/l

Component Information

| Chemical name | Oral LD50 | Dermal LD50 | Inhalation LC50 |
|-----------------------|--------------------|--------------------|-----------------------|
| N,N-Dimethylformamide | = 2800 mg/kg (Rat) | = 1100 mg/kg (Rat) | > 5.85 mg/L (Rat) 4 h |
| 68-12-2 | | | |

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation May cause skin irritation.

Serious eye damage/eye irritation Classification based on data available for ingredients. Causes serious eye irritation.

Respiratory or skin sensitization No information available.

Germ cell mutagenicity No information available.

Carcinogenicity Contains a known or suspected carcinogen. Classification based on data available for

ingredients. May cause cancer.

The table below indicates whether each agency has listed any ingredient as a carcinogen.

| The table below indicates whether each agency has listed any ingredient as a carollogen. | | | | |
|--|-------|----------|-----|------|
| Chemical name | ACGIH | IARC | NTP | OSHA |
| N,N-Dimethylformamide | A3 | Group 2A | - | X |
| 68-12-2 | | | | |

Legend

ACGIH (American Conference of Governmental Industrial Hygienists)

A3 - Animal Carcinogen

IARC (International Agency for Research on Cancer)

Group 2A - Probably Carcinogenic to Humans

OSHA (Occupational Safety and Health Administration of the US Department of Labor)

X - Present

Reproductive toxicity Contains a known or suspected reproductive toxin. Classification based on data available

for ingredients. May damage fertility or the unborn child.

STOT - single exposure No information available.

STOT - repeated exposure No information available.

Aspiration hazard No information available.

Other adverse effects No information available.

Interactive effects No information available.

12. Ecological information

Ecotoxicity

| Chemical name | Algae/aquatic plants | Fish | Toxicity to | Crustacea |
|-----------------------|----------------------|------------------------|----------------|------------------------|
| | | | microorganisms | |
| N,N-Dimethylformamide | EC50: >500mg/L (96h, | LC50: =6300mg/L (96h, | - | EC50: =7500mg/L (48h, |
| 68-12-2 | Desmodesmus | Lepomis macrochirus) | | Daphnia magna) |
| | subspicatus) | LC50: =9800mg/L (96h, | | EC50: =8485mg/L (48h, |
| | | Oncorhynchus mykiss) | | Daphnia magna) |
| | | LC50: =10410mg/L (96h, | | EC50: 6800 - 13900mg/L |
| | | Pimephales promelas) | | (48h, Daphnia magna) |

Persistence and degradability No information available.

Bioaccumulation There is no data for this product.

Component Information

| Chemical name | | Partition coefficient | |
|---------------|-----------------------|-----------------------|--|
| | N,N-Dimethylformamide | -1.028 | |
| | 68-12-2 | | |

Other adverse effects No information available.

13. Disposal considerations

Waste treatment methods

Waste from residues/unused

products

Dispose of in accordance with local regulations. Dispose of waste in accordance with

environmental legislation.

Contaminated packaging Do not reuse empty containers.

14. Transport information

DOT Not regulated

15. Regulatory information

International Inventories

TSCA -.

*Contact supplier for details. One or more substances in this product are either not listed on the US TSCA inventory, listed on the confidential US TSCA inventory or are otherwise exempted from inventory listing requirements

DSL/NDSL -. EINECS/ELINCS -. ENCS -. ECSC -. KECL -. PICCS -. AICS -.

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

| Chemical name | SARA 313 - Threshold Values % |
|---------------------------------|-------------------------------|
| N,N-Dimethylformamide - 68-12-2 | 0.1 |

SARA 311/312 Hazard Categories

Should this product meet EPCRA 311/312 Tier reporting criteria at 40 CFR 370, refer to Section 2 of this SDS for appropriate classifications.

CWA (Clean Water Act)

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302).

| Chemical name | Hazardous Substances RQs | Extremely Hazardous Substances RQs | Reportable Quantity (RQ) |
|----------------------------------|--------------------------|------------------------------------|---|
| N,N-Dimethylformamide 68-12-2 | 100 lb | - | RQ 100 lb final RQ RQ 45.4 kg final RQ |

US State Regulations

California Proposition 65

This product contains the following Proposition 65 chemicals:.

| Chemical name | California Proposition 65 |
|---------------------------------|---------------------------|
| N,N-Dimethylformamide - 68-12-2 | Carcinogen |

U.S. State Right-to-Know Regulations

| Chemical name | New Jersey | Massachusetts | Pennsylvania |
|-----------------------|------------|---------------|--------------|
| N,N-Dimethylformamide | X | X | X |
| 68-12-2 | | | |

U.S. EPA Label Information

EPA Pesticide Registration Number Not applicable

16. Other information

NFPA Health hazards 2 Instability 0 Special hazards -Flammability 1 Health hazards 2 * **HMIS** Flammability 1 Physical hazards 0 Personal protection X

Chronic Hazard Star Legend

* = Chronic Health Hazard

Key or legend to abbreviations and acronyms used in the safety data sheet

Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

TWA Time weighted average Short term exposure limit Ceiling Maximum limit value Skin designation

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR) U.S. Environmental Protection Agency ChemView Database

European Food Safety Authority (EFSA)

EPA (Environmental Protection Agency)

Acute Exposure Guideline Level(s) (AEGL(s))

U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act

U.S. Environmental Protection Agency High Production Volume Chemicals

Food Research Journal

Hazardous Substance Database

International Uniform Chemical Information Database (IUCLID)

Japan GHS Classification

Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)

National Library of Medicine's ChemID Plus (NLM CIP)

National Library of Medicine's PubMed database (NLM PUBMED)

U.S. National Toxicology Program (NTP)

New Zealand's Chemical Classification and Information Database (CCID)

Organization for Economic Co-operation and Development Environment, Health, and Safety Publications

Organization for Economic Co-operation and Development High Production Volume Chemicals Program

Organization for Economic Co-operation and Development Screening Information Data Set

World Health Organization

Revision Date 2023-12-28

Revision NoteNo information available.

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet