

Lenti-X™ Provirus Quantitation Components

Catalog No.	Amount	Lot Number
631240 (Not sold separately) Sold as part of 631239	Each	Specified on product label.

Description

This package contains reagents to be used with the Lenti-X Provirus Quantitation Kit (Cat. No. 631239).

Package Contents

- 30 µl Lenti-X Provirus Control Template (5×10^5 copies/µl)
- 100 µl Lenti-X Provirus Forward Primer (10 µM)
- 100 µl Lenti-X Provirus Reverse Primer (10 µM)
- 4 tubes EASY Dilution Buffer (1 ml/tube)

Storage Conditions

- Store at -20°C .

Expiration Date

- Specified on product label.

Shipping Conditions

- Dry ice

Product Documents

Documents for our products are available for download at takarabio.com/manuals
The following documents apply to this product:

- Lenti-X Provirus Quantitation Kit User Manual

Quality Control Data

Dilutions of the Provirus Control Template and a no-template control (NTC) were tested using the components of the Lenti-X Provirus Quantitation Kit. The resulting standard curve yielded an $R^2 > 0.99$, with an amplification efficiency between 0.9–1.1.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada
800.662.2566
(022423)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999

Lenti-X™ Provirus Quantitation Components

CATALOG NO.

631240

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

2/24/2023