In the interest of conserving resources, we are no longer shipping manuals with products. Please visit www.clontech.com/manuals to obtain an electronic version.

PRODUCT: Human Brain, corpus callosum Poly A* RNA

CATALOG No.

AMOUNT

636133 5 µg

LOT NUMBER 6100206

CONCENTRATION

1 mg/ml in DEPC-treated water

STORAGE CONDITIONS

Store at -70°C

SHELF LIFE

1 year from date of receipt under proper storage conditions.

SHIPPING CONDITIONS

Dry ice (-70°C)

DESCRIPTION

RNA isolated by a modified guanidinium thiocyanate method (1) followed by poly A+ RNA selection with two rounds of oligo(dT)-cellulose columns (2).

POLY A+ RNA SOURCE

normal brain (corpus callosum) pooled from 47 male/female Caucasians, ages 19-76; cause of death: trauma

IMPORTANT NOTE

To prevent contamination by RNases, always wear gloves when handling RNA. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY

QUALITY CONTROL DATA

After electrophoresis on a denaturing gel, the poly A⁺ RNA size ranged from 0.2 – 10 kb.



United States/Canada 800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880 Japan

+81.(0)77.543.6116

Clontech Laboratories, Inc. A Takara Bio Company 1290 Terra Bella Ave. Mountain View, CA 94043 Technical Support (US) E-mail: tech@clontech.com www.clontech.com

REFERENCES

- 1. Ausubel, F. M., et al. (1994) Current Protocols in Molecular Biology. (John Wiley & Sons, NY), Section 4.
- 2. Sambrook, J., et al. (1989) Molecular Cloning: A Laboratory Manual, Second Edition (Cold Spring Harbor Laboratory, Cold Spring Harbor, NY), pp. 6.22-6.34.

APPROVED BY:



Notice to Purchaser

Clontech products are to be used for research purposes only. They may not to be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2006

(PA07059)