

Certificate of analysis

Product Name/Description: TaKaRa MiniBEST Agarose Gel DNA Extraction Kit Ver.4.0

Cat. #: 9762

Lot #: Specified on product label

Storage Conditions: Room temperature

Expiration Date: Specified on product label or package

Package Size: 50 Rxns

Package Contents:	1. Buffer GM	50 ml	
	2. Buffer WB	24 ml	
	3. Elution Buffer	2 ml	x 2
	4. Spin Column		x 50
	5. Collection Tube		x 50

Quality Control Data:

The recovery of DNA fragment (100 bp or 2 kbp) from the DNA band in agarose gel by using this product according to the user manual was more than 50%.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department of Takara Bio Inc.

Notice To Purchaser:

Our products are to be used for Research Use Only. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Product Documents, Patent Markings, Trademarks, & Safety Information:

Referenced below, or available to download from the website or by contacting your local office:
www.takara-bio.co.jp (Japanese)
www.takarabio.com (English)

Licence Statements:

Please visit our website and click the license statement number listed below located on the license statement page.

©2019 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.