

Certificate of Analysis



10X Titanium® *Taq* PCR Buffer

Catalog No(s).

639142

Amount

2,000 rxns

Lot Number

1712366A

Description

10X PCR Buffer for use with the Titanium *Taq* PCR Kit (Cat No. 639210) or Titanium *Taq* DNA Polymerase (Cat. Nos. 639208 & 639209). Buffer is supplied at 5 µl for each 50-µl PCR reaction.

Package Contents

- 10 ml 10X Titanium *Taq* PCR Buffer

Storage Conditions

- Store at –20°C.

Expiration Date

- 01/05/2020

Date of Manufacture

- 01/05/2018

Shipping Conditions

- Dry ice (–70°C)

Product Documents

Documents for our products are available for download at takarabio.com/manuals

Quality Control Data

The 10X Titanium *Taq* PCR Buffer was tested using Titanium *Taq* DNA Polymerase, and other reagents from the Titanium *Taq* PCR Kit, in a 50-µl PCR reaction using 5 µl of Marathon®-Ready Human Placenta cDNA (Cat. No. 639311) and control primers specific for a 1.3-kb fragment from the transferrin receptor gene (0.2 µM each). Conditions were set at:

1 cycle	95°C, 1 min
30 cycles	95°C, 15 sec
	68°C, 1.5 min

5 µl of the PCR product was electrophoresed on a 1.2% TAE/agarose gel to confirm the presence of a 1.3-kb band with minimal background.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.comUnited States/Canada
800.662.2566
(111616)Asia Pacific
+1.650.919.7300Europe
+33.(0)1.3904.6880Japan
+81.(0)77.565.6999

10X Titanium® Taq PCR Buffer

CATALOG NO.

639142

NOTICE TO PURCHASER:

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

1/25/2018