

Certificate of Analysis

DO Supplement –His/–Ura

Catalog No.

630422

Amount

10 g

Lot Number

Specified on product label.

Description

Dropout (DO) supplement designed to be used at a concentration of 0.76 g/L with our Minimal SD Base (Cat. No. 630411) or Minimal SD Agar Base (Cat. No. 630412) to create synthetic, defined medium lacking the specified nutrients. Please see our Yeast Media Protocol-At-A-Glance for details.

Package Contents

- 10 g DO Supplement –His/–Ura

Storage Conditions

- Store at room temperature

Shelf Life

- Specified on product label

Shipping Conditions

- Room temperature

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Yeast Media Protocol-At-A-Glance

Quality Control Data

The following inspections were completed and confirmed:

- Powder appearance - free-flowing, very fine/smooth
- Powder color - white

The powder was dissolved at a concentration of 0.76 g/L. The appearance, color, and pH were as expected.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.comUnited States/Canada
800.662.2566
(090121)Asia Pacific
+1.650.919.7300Europe
+33.(0)1.3904.6880Japan
+81.(0)77.565.6999

DO Supplement –His/–Ura

CATALOG NO.

630422

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

9/1/2021