

Magnetic DYKDDDDK Immunoprecipitation Kit

Catalog No.	Amount
635694	50 rxns

Description

The Magnetic DYKDDDDK Immunoprecipitation Kit contains the monoclonal anti-DYKDDDDK antibody, which recognizes the DYKDDDDK tag, crosslinked to magnetic beads. This product may be used to purify or immunoprecipitate DYKDDDDK-tagged fusion proteins from cell lysates. The amount of crosslinked antibody on the magnetic beads is 8 mg/ml. The beads are provided as 25% slurry in storage buffer. The kit also includes the Magnetic Beads Immunoprecipitation Buffer Set and the ProteoGuard™ EDTA-Free Protease Inhibitor Cocktail.

Package Contents

- Anti-DYKDDDDK Magnetic Beads (Cat. No. 635695) >> [View Components](#)
- Magnetic Beads Immunoprecipitation Buffer Set (Cat. No. 635696) >> [View Components](#)
- ProteoGuard™ EDTA-Free Protease Inhibitor Cocktail (Cat. No. 635673) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for Clontech® products are available for download at www.clontech.com/manuals
The following documents apply to this product:

- Magnetic DYKDDDDK Immunoprecipitation Kit Protocol-At-A-Glance

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech, the Clontech logo, and ProteoGuard are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2013 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (091313)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116