

PRODUCT COMPONENTS LIST: Yeast Media Set 1 Plus



CATALOG No. 630493

DESCRIPTION

The Yeast Media Set 1 Plus contains the media and supplements needed for yeast one-hybrid library screening with the Matchmaker™ Gold Yeast One-Hybrid Library Screening System (Cat. No. 630491). For the ready-to-go, preformulated pouches; just add water and autoclave according to the instructions provided in the Clontech Yeast Media Protocol-at-a-Glance (PT4057-2). No measuring, mixing, or pH adjustments are required. Directions for the use of Aureobasidin A are provided in the Matchmaker Gold Yeast One-Hybrid Library Screening System User Manual (PT4087-1).

PACKAGE CONTENTS

- **1 each Yeast Media Set 1 (Cat. No. 630492)**
 - 2 pouches YPDA Broth
 - 1 pouch YPDA with Agar
 - 1 pouch SD/-Leu Broth
 - 10 pouches SD/-Leu with Agar
 - 2 pouches SD/-Ura with Agar
- **1 mg Aureobasidin A (Cat. No. 630466)**

OTHER

- Clontech Yeast Media Protocol-at-a-Glance (PT4057-2)

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

For licensing statements and storage conditions, please see the Certificate of Analysis supplied with each component.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2008 Clontech Laboratories, Inc.

(PA8Y2981)