

In the interest of conserving resources, we are no longer shipping manuals with products. Please visit [www.clontech.com/manuals](http://www.clontech.com/manuals) to obtain an electronic version.

## PRODUCT COMPONENTS LIST: Lenti-X™ Maxi Purification Kit (with Rack)

CATALOG No. 631245

### DESCRIPTION

The Lenti-X Maxi Purification Kit purifies lentivirus from crude supernatants, and produces outstanding yields of highly purified virus. The gravity column-based protocol is fast, simple, and effective. The Lenti-X column gently retains charged virus particles from the supernatant, while unbound material exits the column during the wash steps. The purified virus that is eluted from the column is recovered fully intact and fully functional. The convenient Collection Rack L holds up to 3 Lenti-X columns and 50 ml conical tubes, and secures them in place for spill-free purifications.

### PACKAGE CONTENTS

- 1 each Lenti-X Maxi Purification Kit (5 preps) (Cat. No. 631234)
- 1 each Collection Rack L (Cat. No. 631246)

### USER MANUALS

User manuals for Clontech products are available for download at [www.clontech.com/manuals](http://www.clontech.com/manuals). The following user manuals apply to this product:

- Lenti-X Maxi Purification Kit Protocol-at-a-Glance (PT5024-2)



**Clontech**

**United States/Canada**  
800.662.2566

**Asia Pacific**  
+1.650.919.7300

**Europe**  
+33.(0)1.3904.6880

**Japan**  
+81.(0)77.543.6116

Clontech Laboratories, Inc.  
A Takara Bio Company  
1290 Terra Bella Ave.  
Mountain View, CA 94043  
Technical Support (US)  
E-mail: [tech@clontech.com](mailto:tech@clontech.com)  
[www.clontech.com](http://www.clontech.com)

### Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

For licensing statements and storage conditions, please see the Certificate of Analysis supplied with each component.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.

(PA9Z3498)