

Guide-it™ Cas9 Monoclonal Antibody (Clone TG8C1)

Catalog No.	Amount
632627	3 x 100 µg

Description

A monoclonal antibody produced by hybridoma cells against full-length recombinant Cas9 protein from *Streptococcus pyogenes*. This antibody recognizes as little as 1 ng of wild-type Cas9 in bacterial and mammalian cell lysates by Western blot.

Package Contents

- 3 x Guide-it Cas9 Monoclonal Antibody (Clone TG8C1) (Cat. No. 632628) >>[View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for Clontech® products are available for download at www.clontech.com/manuals

The following documents apply to this product:

- Antibodies Protocol Guide (PT3407-1)

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech®, the Clontech logo, and Guide-it are trademarks of Clontech Laboratories, Inc. All other trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech Laboratories, Inc. is a Takara Bio Company. ©2016 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (072116)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------