

SMARTer® Stranded Total RNA Sample Prep Kit - Low Input Mammalian

Catalog No.	Amount
634861	24 rxns

Description

The SMARTer Stranded Total RNA Sample Prep Kit - Low Input Mammalian provides a complete solution for sensitive and streamlined total RNA analysis. The kit is designed to work with a wide range of input amounts (10–100 ng) of total RNA of any quality. A variety of human, rat, and mouse samples can be used, including total RNA from FFPE material, tissue samples, etc. The kit provides whole transcriptome analysis with accurate measurement of strand orientation, unbiased coverage, high transcript mapping, and gene counts.

Package Contents

- SMARTer Stranded RNA-Seq Components (Cat. No. 634841; Not sold separately) >>[View Components](#)
- SeqAmp™ DNA Polymerase (Cat. No. 638504) >>[View Components](#)
- 2 x Illumina® Indexing Primer Set (Cat. No. 634844; Not sold separately) >>[View Components](#)
- RiboGone™ - Mammalian (Cat. No. 634847) >>[View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for Clontech® products are available for download at www.clontech.com/manuals

The following documents apply to this product:

- RiboGone - Mammalian Protocol-At-A-Glance
- SMARTer Stranded RNA-Seq Kit User Manual
- SeqAmp DNA Polymerase Protocol-At-A-Glance

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Illumina is a trademark of Illumina, Inc.

Clontech, the Clontech logo, RiboGone, SeqAmp, and SMARTer are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2014 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (021014)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.543.6116
--	---------------------------------	------------------------------	-----------------------------