

## HisTALON™ Superflow Cartridge Purification Kit

Catalog No.	Amount
635649	5 x 1 ml cartridges

### Description

Prepacked cartridges, each containing 1 ml of TALON® Superflow Resin, an immobilized metal affinity resin for the purification of recombinant his-tagged proteins. These cartridges may be used under native or denaturing conditions, in combination with a set of equilibration and elution buffers optimized for purifying his-tagged proteins using HisTALON Cartridges.

### Package Contents

- HisTALON Superflow Cartridges (5 x 1 ml cartridges)  
(Cat. No. 635650) >> [View Components](#)
- HisTALON Buffer Set (20 purifications using 1 ml cartridges)  
(Cat. No. 635651) >> [View Components](#)

### Product Documents

Documents for Clontech products are available for download at [www.clontech.com/manuals](http://www.clontech.com/manuals).

The following document applies to this product:

- HisTALON Superflow Cartridges User Manual (PT4094-1)

---

### Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech, the Clontech logo, HisTALON, and TALON are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2012 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

---

### Clontech Laboratories, Inc.

A Takara Bio Company  
1290 Terra Bella Avenue, Mountain View, CA 94043, USA  
U.S. Technical Support: [tech@clontech.com](mailto:tech@clontech.com)

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (030112)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116