

**PRODUCT COMPONENTS LIST: His60 Ni Gravity Column Purification Kit****CATALOG No.** 635658**DESCRIPTION**

The His60 Ni Gravity Column Purification Kit provides prepacked, high capacity His60 Ni Gravity Columns and all of the buffers necessary to purify recombinant his-tagged proteins. This kit is suitable for gravity-flow purification under native or denaturing conditions. Each column contains 1 ml of His60 Ni Superflow™ Resin, and can be reused 4 times to purify different preparations of the same protein. The kit provides enough reagents for 20 purifications.

**PACKAGE CONTENTS**

- His60 Ni Gravity Columns (5 columns)  
(Cat. No. 635657)
- His60 Ni Buffer Set  
(Cat. No. 635665)

*Note: For more information, please see the Certificate of Analysis supplied with each component.*

**Clontech****United States/Canada**  
800.662.2566**Asia Pacific**  
+1.650.919.7300**Europe**  
+33.(0)1.3904.6880**Japan**  
+81.(0)77.543.6116

Clontech Laboratories, Inc.  
A Takara Bio Company  
1290 Terra Bella Ave.  
Mountain View, CA 94043  
Technical Support (US)  
E-mail: tech@clontech.com  
www.clontech.com

**Notice to Purchaser**

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Superflow™ is a trademark of Sterogene Bioseparations, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2009 Clontech Laboratories, Inc

(PA933158)