

# Product Components List

## In-Fusion® HD Cloning Plus CE

Catalog No(s).	Amount
638916	10 rxns
638917	50 rxns
638918	100 rxns

### Description

The In-Fusion HD Cloning Plus CE kit enables high-efficiency, high-fidelity, directional cloning of one or more PCR fragments into any vector. In addition to the In-Fusion HD Cloning Kit, this package includes: Cloning Enhancer, which removes background plasmid DNA and PCR residue, eliminating the need for PCR insert purification prior to cloning; Stellar™ Competent Cells, an *E. coli* HST08 strain that provides high transformation efficiency paired with blue-white screening capability when used with pUC plasmid vectors; and CloneAmp™ HiFi PCR Premix, a convenient high-fidelity PCR 2X master mix that provides exceptionally accurate and efficient DNA amplification.

### Package Contents

#### Cat. No. 638916 (10 rxns):

- In-Fusion HD Cloning Kit (Cat. No. 639648) >> [View Components](#)
- Cloning Enhancer (Cat. No. 639613) >> [View Components](#)
- Stellar Competent Cells (Cat. No. 636763) >> [View Components](#)
- CloneAmp HiFi PCR Premix (Cat. No. 638500) (Not sold separately) >> [View Components](#)

#### Cat. No. 638917 (50 rxns):

- In-Fusion HD Cloning Kit (Cat. No. 639649) >> [View Components](#)
- Cloning Enhancer (Cat. No. 639614) >> [View Components](#)
- 5 x Stellar Competent Cells (Cat. No. 636763) >> [View Components](#)
- 2 x CloneAmp HiFi PCR Premix (Cat. No. 638500) (Not sold separately) >> [View Components](#)

#### Cat. No. 638918 (100 rxns):

- In-Fusion HD Cloning Kit (Cat. No. 639650) >> [View Components](#)
- Cloning Enhancer (Cat. No. 639615) >> [View Components](#)
- 10 x Stellar Competent Cells (Cat. No. 636763) >> [View Components](#)
- CloneAmp HiFi PCR Premix (Cat. No. 639298) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

---

#### Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA  
U.S. Technical Support: [techUS@takarabio.com](mailto:techUS@takarabio.com)

United States/Canada 800.662.2566 (102318)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------

# Product Components List

Cat. No(s). 638916, 638917 & 638918

In-Fusion® HD Cloning Plus CE

---

## Product Documents

Documents for our products are available for download at [takarabio.com/manuals](http://takarabio.com/manuals)

The following documents apply to this product:

- In-Fusion HD Cloning Kit User Manual (PT5162-1)
  - pUC19 Linearized Vector Information (PT4065-5)
  - Stellar Competent Cells Protocol-at-a-Glance (PT5055-2)
  - CloneAmp HiFi PCR Premix Protocol-At-A-Glance
- 

## Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at [takarabio.com](http://takarabio.com). It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at [takarabio.com](http://takarabio.com).

This document has been reviewed and approved by the Quality Department.