

Product Components List

SMART-Seq® Single Cell PLUS Kit

Catalog Nos.	Amount
R400750	48 rxns
R400751	96 rxns

Description

The SMART-Seq Single Cell PLUS kit is a complete kit designed to generate high-quality, Illumina® sequencing-ready libraries from single cells known to have low RNA content. The SMART-Seq Single Cell PLUS Kit protocol involves generating cDNA using the SMART-Seq Single Cell Kit, fragmenting the cDNA using the SMART-Seq Library Prep Kit (not sold separately), and adding indexes using the unique dual index kit.

Package Contents

R400750 (48 rxns):

- SMART-Seq Single Cell Kit (Cat. No. 634471) >>[View component](#)
- SMART-Seq Library Prep Kit (not sold separately; Cat. No. R400746) >>[View component](#)
- Unique Dual Index Kit - 48U (Cat. No. R400744) >>[View component](#)

R400751 (96 rxns):

- SMART-Seq Single Cell Kit (Cat. No. 634472) >>[View component](#)
- SMART-Seq Library Prep Kit (not sold separately; Cat. No. R400747) >>[View Component](#)
- Unique Dual Index Kit - 96U (Cat. No. R400745) >>[View Component](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- SMART-Seq Single Cell PLUS Kit User Manual
- Unique Dual Index Kit Protocol-At-A-Glance

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2020 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: techUS@takarabio.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (030520)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999