

pGBKm Vector

Catalog No(s).

631851




Figure 1. pGBKm vector map.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (062113)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.543.6116
--	---------------------------------	------------------------------	-----------------------------

Vector Map

pGBKm

631851

Location of Features

- TRP1 (isomerase gene required for tryptophan biosynthesis): 103–777 (complementary)
- pTRP1 (TRP1 promoter): 778–1059 (complementary)
- 2 μ ori (yeast 2 μ origin of replication): 1318–2660
- pAmpR Promoter (ampicillin resistance promoter): 2687–2791
- NeoR/KanR (Tn5 neomycin/kanamycin resistance gene): 2792–3586
- pUC ori (pUC origin of replication): 4194–4782

Notice to Purchaser

Clontech® products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

Clontech and the Clontech logo are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2013 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.