

Lenti-X™ Accelerator Starter Kit

Catalog No(s).	Amount
631254	1 kit

Description

Lenti-X Accelerator is a magnetic bead-based technology designed to accelerate lentiviral transduction experiments. Traditional lentiviral transductions require overnight incubation of your cells and virus with Polybrene®, but Lenti-X Accelerator allows you to complete this process in only 30 minutes and avoid using Polybrene altogether. Charged magnetic beads bind to virus particles from your lentiviral packaging supernatant and are pulled into contact with your cells in the presence of a magnetic field. The Lenti-X Accelerator reagent and a Magnetic Separator for Cell Culture are included in the Lenti-X Accelerator Starter Kit. Lenti-X Accelerator is excellent for transducing lentivirus, MMLV retrovirus, and MSCV retrovirus preparations.

Package Contents

- 400 µl Lenti-X Accelerator (Cat. No. 631256) >>[View Components](#)
- 1 each Magnetic Separator for Cell Culture (Cat. No. 631255)

Product User Manuals

User manuals for Clontech products are available for download at www.clontech.com/manuals

The following user manuals apply to this product:

- Lenti-X Accelerator Protocol (PT5159-2)

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Polybrene® is a registered trademark of Abbott Laboratories Corp.

Clontech®, the Clontech logo, and Lenti-X are trademarks of Clontech Laboratories, Inc. All other trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech Laboratories, Inc. is a Takara Bio Company. ©2016 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (053116)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999