

Product Components List

BacPAK™ qPCR Titration Kit

Catalog No.	Amount
631414	1 kit

Description

This kit can be used along with real-time PCR equipment to rapidly quantify the virion content of baculoviral supernatants. The copy numbers of viral DNA genomes contained in virus samples are determined by comparing the qPCR threshold cycles (C_t) derived from serial dilutions of purified viral DNA to the C_t values derived from dilutions of a calibrated control template. The kit includes a viral nucleic acid purification kit, a calibrated BacPAK DNA Control Template, primers, and all the reagents required to conduct 200 qPCR reactions using either the BacPAK DNA Control Template or purified viral DNA samples.

Package Contents

- BacPAK qPCR Titration Kit Components (Cat. No. 631415) (Not sold separately) >> [View Components](#)
- TB Green® Advantage® qPCR Premix (200 rxns) (Cat. No. 639676) >> [View Components](#)
- NucleoSpin Virus (Cat. No. 740983.10) (Not sold separately) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals
The following documents apply to this product:

- BacPAK qPCR Titration Kit User Manual

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: techUS@takarabio.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (050619)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999