

Product Components List

Guide-it™ Indel Identification Kit

Catalog Nos.	Amount
631444	10 Rxns

Description

The Guide-it Indel Identification Kit is used for characterization of insertions and deletions (indels) generated by gene editing tools, such as CRISPR/Cas9. This kit contains all components needed to amplify, clone, and prepare modified target sites for DNA sequence analysis. This kit uses Terra™ PCR Direct to amplify targets directly from crude genomic DNA extracts. The resulting pool of fragments, which may contain a variety of indels, is cloned into a prelinearized pUC19 vector using the In-Fusion® cloning system. Colony PCR of individual clones using Terra PCR Direct followed by DNA sequencing allows indel characterization.

Package Contents

- Guide-it Indel Identification Components (Cat. No. 631445) (Not sold separately) >> [View Components](#)
- In-Fusion Snap Assembly Master Mix (Cat. No. 638947) >> [View Components](#)
- NucleoSpin Gel and PCR Clean-Up Kit (Cat. No. 740609.10) >> [View Components](#)
- Stellar™ Competent Cells (Cat. No. 636763) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Guide-it Indel Identification Kit User Manual
- NucleoSpin PCR Clean-Up & Gel Extraction User Manual
- pUC19 Linearized Vector Information

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2021 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: technical_support@takarabio.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (042121)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999