

In the interest of conserving resources, we are no longer shipping manuals with products. Please visit www.clontech.com/manuals to obtain an electronic version.

PRODUCT COMPONENTS LIST: Pluripotency Check PCR Kit

CATALOG No.	SIZE
631965	each

DESCRIPTION

The Pluripotency Check PCR Kit can be used to validate the pluripotent status (ability to differentiate into specialized cells of all three germ layers) of your stem cell culture. It includes 16 forward and reverse primer mixes, based on nine genes expressed in pluripotent stem cells and two control genes, plus the reagents necessary for RNA extraction, first-strand cDNA synthesis, and PCR.

PACKAGE CONTENTS

- Pluripotency Check PCR Primer Set (16 x 25 rxns)
(Cat. No. 631966)
- NucleoSpin[®] RNA II (10 preps)
(Cat. No. 740955.10)
- SMART[™] MMLV Reverse Transcriptase (8,000 units)
(Cat. No. 639523)
- TITANIUM[®] Taq PCR Kit (30 rxns)
(Cat. No. 639211)

Note: for more information, please see the Certificate of Analysis supplied with each component.



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2009 Clontech Laboratories, Inc.

(PA993414)