

Product Components List

DD-ZsGreen1 Reporter System

Catalog No.	Amount
632192	each

Description

The DD-ZsGreen1 Reporter System includes the pDD-ZsGreen1 Reporter vector and Shield1. The pDD-ZsGreen1 Reporter is a promoterless vector designed to insert your promoter of interest upstream of the green fluorescent protein ZsGreen1, tagged at its N-terminus with the ProteoTuner™ destabilization domain (DD). The DD causes the reporter protein to be rapidly targeted to and degraded by proteasomes. This very efficient and controllable destabilization method will minimize any reporter background fluorescence from leaky promoters prior to promoter activation.

To analyze promoter activity, a candidate inducer is added to the medium along with Shield1, which effectively stabilizes the reporter protein and allows it to accumulate. As a result, only the reporter molecules expressed during promoter induction will contribute to the fluorescence signal, providing a considerably higher signal-to-noise ratio than that obtained with non-destabilized or constitutively destabilized reporter systems.

Package Contents

- pDD-ZsGreen1 Reporter (Cat. No. 632195) (Not sold separately) >> [View Components](#)
- Shield1 (500 µl) (Cat. No. 632189) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at www.clontech.com/manuals

The following documents apply to this product:

- ProteoTuner Plasmid-Based Shield Systems User Manual
- DD-Fluorescent Protein Reporter Systems Protocol-at-a-Glance (PT4088-2)

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

©2016 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (102616)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------