

Product Components List

Guide-it™ Complete sgRNA Screening System

Catalog No.	Amount
632636	50 rxns

Description

The Guide-it Complete sgRNA Screening System is used to both synthesize and test the efficacy of single guide RNAs (sgRNAs) for studies involving Cas9-mediated gene editing. The kit contains all the components needed to both *in vitro* transcribe sgRNAs and test their efficacy against a DNA target fragment using recombinant Cas9 nuclease.

Package Contents

- Guide-it sgRNA *In Vitro* Transcription Components v2 (Cat. No. 632637; Not sold separately) [>>View Components](#)
- Guide-it IVT RNA Clean-Up Kit (Cat. No. 632638) [>>View Components](#)
- Guide-it sgRNA Screening Kit (Cat. No. 632639) [>>View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at www.clontech.com/manuals

The following documents apply to this product:

- Guide-it sgRNA *In Vitro* Transcription and Screening Systems User Manual

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

©2016 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at <http://www.takarabio.com>.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (102016)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------