

# Product Components List

## Embgenix™ PGT-A Kit (RUO)

**Catalog No.**  
634760

**Amount**  
96 rxns

### Description

Embgenix PGT-A Kit (RUO) has been developed to analyze DNA obtained from biopsies of single or multiple embryonic cells for preimplantation genetic testing for aneuploidies (PGT-A). All 24 chromosomes are accurately screened for whole chromosome aneuploidy and sub-chromosomal abnormalities utilizing whole genome amplification (WGA) and next generation sequencing (NGS) technologies.

### Package Contents

- Embgenix PGT-A Core Kit (RUO) (Cat. No. 634761) (Not sold separately) >> [View Components](#)
- Unique Dual Index Kit (1–96) (Cat. No. 634752) >> [View Components](#)
- NucleoMag NGS Clean-up and Size Select (Cat. No. 744970.5) >> [View Product Page](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

### Product Documents

Documents for our products are available for download at [takarabio.com/manuals](http://takarabio.com/manuals)

The following documents apply to this product:

- Embgenix PGT-A Kit (RUO) User Manual for Illumina® MiSeq® System
- Embgenix PGT-A Kit (RUO) Protocol-At-A-Glance for Illumina MiSeq System
- Embgenix PGT-A Kit (RUO) Library Pooling Calculator for Illumina MiSeq System
- Embgenix PGT-A Kit (RUO) User Manual for Illumina NextSeq® 500/550 System
- Embgenix PGT-A Illumina MiSeq96 SampleSheet LocalRunManager
- Embgenix PGT-A Illumina MiSeq96 SampleSheet
- Embgenix PGT-A Illumina NextSeq96 SampleSheet
- Embgenix Analysis Software (RUO) User Manual

### Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at [takarabio.com](http://takarabio.com). It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2022 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at [takarabio.com](http://takarabio.com).

This document has been reviewed and approved by the Quality Department.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

United States/Canada 800.662.2566 (080322)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------