

## Marathon® cDNA Amplification Kit

Catalog No.	Amount
634913	100 rxns

### Description

This kit contains reagents for (1) making adaptor-ligated ds cDNA from poly A+ RNA and (2) performing 5' and 3' Marathon RACE PCR using gene-specific primers and a mixture of thermostable DNA polymerases supplied by the user. This kit contains sufficient reagents for 5 cDNA synthesis reactions and 100 PCR reactions.

### Package Contents

- Marathon cDNA Amplification Kit (Cat. No. 634952) (Not sold separately) [>>View Components](#)
- Nucleotrap Gel Extraction Kit (Cat. No. 740584.10) (Not sold separately) [>>View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

### Product Documents

Documents for Clontech® products are available for download at [www.clontech.com/manuals](http://www.clontech.com/manuals)

The following documents apply to this product:

- Marathon cDNA Amplification Kit User Manual (PT1115-1)
- Marathon cDNA Amplification Kit Protocol-at-a-Glance (PT1115-2)
- NucleoTraPCR PCR Clean-up & NucleoTrap Gel Extraction User Manual

---

### Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech, the Clontech logo, and Marathon are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2015 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

---

### Clontech Laboratories, Inc.

A Takara Bio Company  
1290 Terra Bella Avenue, Mountain View, CA 94043, USA  
U.S. Technical Support: [tech@clontech.com](mailto:tech@clontech.com)

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (051115)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116