

D/D Solubilizer

Catalog No.	Amount
635053	5 x 500 µl

Description

The D/D Solubilizer is a membrane-permeant ligand that solubilizes complexes formed by proteins fused to the self-assembly domain (DmrD), provided in the iDimerize™ Reverse Dimerization System (635066).

Package Contents

- 5 tubes D/D Solubilizer (Cat. No. 635054) [>> View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product User Manuals

User manuals for Clontech products are available for download at www.clontech.com/manuals.

The following user manuals apply to this product:

- iDimerize Reverse Dimerization System User Manual (PT5180-1)

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

These products are sold under license from ARIAD Pharmaceuticals, Inc.

Clontech, the Clontech logo, and iDimerize are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2011 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116
(PA134456)			