

Product Components List

A/C Heterodimerizer

Catalog No.	Amount
635056	5 x 500 µl

Description

A/C Heterodimerizer is a membrane-permeant ligand that causes heterodimerization of proteins fused to two distinct dimerization domains—DmrA and DmrC—provided in the iDimerize™ Inducible Heterodimer System (635067) and the iDimerize Regulated Transcription System (635081). **NOTE:** The A/C Heterodimerizer is identical to the AP21967 ligand previously supplied in the ARGENT Regulated Heterodimerization Kit from ARIAD Pharmaceuticals, Inc.

Package Contents

- 5 tubes A/C Heterodimerizer (Cat. No. 635057) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- iDimerize Inducible Heterodimer System User Manual (PT5179-1)
- iDimerize Regulated Transcription System User Manual

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

©2017 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (031617)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999