

## iDimerize™ Regulated Transcription System

### Catalog No.

635081

### Description

The iDimerize Regulated Transcription System allows tightly regulated, inducible expression of your gene of interest by the addition of a small molecule to the cell medium. The system includes the iDimerize Regulated Transcription Vector Set 2, which contains a vector that encodes two transcription factor domains (which combine to create the inducible transcription factor), a mammalian expression vector used to express your gene of interest, and a control vector. The system also includes the A/C Heterodimerizer, the membrane-permeant compound required for induction of gene expression through the heterodimerization of the two transcription factor domains.

### Package Contents

- iDimerize Regulated Transcription Vector Set 2 (Cat. No. 635082; Not sold separately) >> [View Components](#)
- 500 µl A/C Heterodimerizer (0.5 mM; Cat. No. 635057) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

### Product Documents

Documents for Clontech products are available for download at [www.clontech.com/manuals](http://www.clontech.com/manuals)

The following documents apply to this product:

- iDimerize Regulated Transcription System User Manual

---

### Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

iDimerize products are sold under license from ARIAD Pharmaceuticals, Inc.

Clontech, the Clontech logo, and iDimerize are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2012 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

---

### Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: [tech@clontech.com](mailto:tech@clontech.com)

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (071812)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116