

A/C Heterodimerizer

Catalog No.	Amount
635095	5 x 5 mg

Description

A/C Heterodimerizer is a membrane-permeant ligand that causes heterodimerization of proteins fused to two distinct dimerization domains—DmrA and DmrC—provided in the iDimerize™ Heterodimer Systems.

NOTE: The A/C Heterodimerizer is identical to the AP21967 ligand previously supplied in the ARGENT Regulated Heterodimerization Kit from ARIAD Pharmaceuticals, Inc.

Package Contents

- 5 x 5 mg A/C Heterodimerizer (Cat. No. 635055) >>[View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for Clontech® products are available for download at www.clontech.com/manuals

The following documents apply to this product:

- iDimerize Inducible Heterodimer System User Manual (PT5179-1)
- iDimerize Inducible Heterodimer System (with Tet-On 3G) User Manual
- iDimerize Regulated Transcription System User Manual
- Lenti-X iDimerize Inducible Heterodimer System User Manual

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

These products are sold under license from ARIAD Pharmaceuticals, Inc.

Clontech®, the Clontech logo, and iDimerize are trademarks of Clontech Laboratories, Inc. All other trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech Laboratories, Inc. is a Takara Bio Company. ©2016 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (052616)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--------------------------------------------------	---------------------------------	------------------------------	-----------------------------