

His60 Ni Superflow Resin & xTractor™ Buffer Bundle

Catalog No(s).	Amount
635676	20 purifications

Description

His60 Ni Superflow Resin is a high-capacity immobilized metal affinity chromatography (IMAC) resin for the efficient purification of recombinant his-tagged proteins under native or denaturing conditions. xTractor Buffer is an optimized buffer system for extraction of proteins from fresh or frozen bacterial cells, mammalian cells, yeast cells, or baculovirus-infected cells. This buffer is useful in extracting proteins over a broad range of molecular weights, including high molecular weight proteins that cannot be extracted unless the membrane is fully disrupted. xTractor Buffer has been optimized for extraction of his-tagged proteins, so it is compatible with all IMAC resins. This buffer is suitable for any scale protein extraction & any tag (use with his-, GST-, FLAG- HA- or other tagged proteins).

Package Contents

- His60 Ni Superflow Resin (25 ml)
(Cat. No. 635660) >>[View Components](#)
- xTractor Buffer (100 ml)
(Cat. No. 635656) >>[View Components](#)

Product Documents

Documents for Clontech® products are available for download at www.clontech.com/manuals.

The following documents apply to this product:

- His60 Ni Superflow Resin & Gravity Column User Manual (PT5017-1)
- xTractor Buffer & xTractor Buffer Kit User Manual (PT3735-1)
- xTractor Buffer Bacterial Sample Preparation Protocol (PT3735-2)

Notice to Purchaser

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech®, the Clontech logo, and xTractor are trademarks of Clontech Laboratories, Inc. All other trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech Laboratories, Inc. is a Takara Bio Company. ©2016 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 (030916)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.543.6116
--	---------------------------------	------------------------------	-----------------------------