

Product Components List

Capturem™ His-Tagged Purification Maxiprep Columns

Catalog No.	Amount
635719	50 columns

Description

Capturem His-Tagged Purification Maxiprep Columns are single-use disposable membrane spin columns for simple, rapid purification of his-tagged proteins in up to 25 ml of clarified lysate. The columns are suitable for use under native or denaturing conditions, and in the presence of additives such as DTT (up to 10 mM), β ME (up to 30 mM), TCEP (up to 5 mM), EDTA (up to 10 mM), or glycerol.

Package Contents

- 2 x Capturem His-Tagged Purification Maxiprep Columns (25 columns; Cat. No. 635715) >>[View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Capturem His-Tagged Purification Maxiprep Columns Protocol-At-A-Glance

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: techUS@takarabio.com

United States/Canada 800.662.2566 (110218)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------