

Product Components List

Prelude™ One-Step PreAmp Master Mix

Catalog No.	Amount
638553	40 rxns

Description

Prelude One-Step PreAmp Master Mix is a 2X RT-PCR mix optimized for unbiased co-amplification of multiple targets (>100) in a single tube starting from low amounts of input RNA (100 pg mRNA or 10 ng total RNA). Amplified products can then be used for downstream qPCR, genotyping, or targeted enrichment analysis.

Package Contents

- 4 x Prelude One-Step PreAmp Master Mix (10 rxns; Cat. No. 638554) >>[View components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Prelude One-Step PreAmp Master Mix Protocol-At-A-Glance

References

Ruano, G., *et al.* Biphasic amplification of very dilute DNA samples via “booster” PCR. *Nucleic Acids Res.* **17**, 5407 (1989).

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product’s web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2020 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada 800.662.2566 (111620)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------