

In the interest of conserving resources, we are no longer shipping manuals with products. Please visit [www.clontech.com/manuals](http://www.clontech.com/manuals) to obtain an electronic version.

## PRODUCT COMPONENTS LIST: Terra™ PCR Direct Polymerase Mix

CATALOG No. 639271

AMOUNT: 800 rxns

### DESCRIPTION

Terra PCR Direct Polymerase Mix is a novel enzyme developed for direct amplification from tissue samples, crude extracts, and dirty templates. It's perfect for amplifying short DNA targets (up to 2 kb), regardless of GC content or template purity. Terra PCR Direct is a highly sensitive enzyme that allows amplification of targets from small amounts of template. Moreover, the enzyme comes pre-blended with a monoclonal antibody that suppresses polymerase activity up to 98°C, allowing automatic hot start PCR.

### PACKAGE CONTENTS

- 4 each Terra PCR Direct Polymerase Mix (200 rxns)  
(Cat. No. 639270)

### PRODUCT USER MANUALS

User manuals for Clontech products are available for download at [www.clontech.com/manuals](http://www.clontech.com/manuals). The following user manuals apply to this product:

- Terra PCR Direct Polymerase Mix User Manual (PT5126-1)
- Terra PCR Direct Polymerase Mix Protocol-at-a-Glance (PT5126-2)

*Note: For more information, please see the Certificate of Analysis supplied with each component.*



**Clontech**

**United States/Canada**  
800.662.2566

**Asia Pacific**  
+1.650.919.7300

**Europe**  
+33.(0)1.3904.6880

**Japan**  
+81.(0)77.543.6116

Clontech Laboratories, Inc.  
A Takara Bio Company  
1290 Terra Bella Ave.  
Mountain View, CA 94043  
Technical Support (US)  
E-mail: [tech@clontech.com](mailto:tech@clontech.com)  
[www.clontech.com](http://www.clontech.com)

### Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.

(PA013513)