

High Fidelity PCR EcoDry™ Premix

Catalog No.	Amount	Lot Number
639280	48 rxns	Specified on product label.

Description

High Fidelity PCR EcoDry Premix is a convenient dry master mix that allows efficient and accurate amplification of DNA. When reconstituted, each tube of High Fidelity PCR EcoDry Premix contains all of the necessary reagents for high fidelity PCR, including Advantage® 2 Polymerase Mix, dNTPs, and optimized buffer. Advantage 2 combines a 5'-exonuclease-deficient *Taq* DNA Polymerase, a hot start antibody, and a small amount of a proofreading polymerase for improved fidelity. Simply reconstitute the mix by adding PCR-grade water along with your primers and template. EcoDry Premixes are packaged as individual 8-well tube strips with optically clear sealing caps.

Package Contents

- 2 High Fidelity PCR EcoDry Premix (24 rxns each; Cat. No. 639282) >> [View Components](#)

Product User Manuals

User manuals for Clontech products are available for download at www.clontech.com/manuals.

The following user manuals apply to this product:

- High Fidelity PCR EcoDry Premix Protocol-at-a-Glance (PT5152-2)

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

For licensing statements and storage conditions, please see the Certificate of Analysis supplied with each component.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (PA063933)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116