

Certificate of Analysis

PrepX™ ILMN Barcodes 1-24 (400077-1)

Catalog No.

640103

Amount

24 barcodes

Lot Number

Specified on product label.

Description

Each plate of PrepX ILMN Barcodes 1-24 allows the preparation of up to 24 Illumina®-ready DNA libraries. Along with a barcode adapter, each well contains a universal adapter for sequencing on Illumina platforms.

Package Contents

Package 1:

- 1 plate PrepX ILMN Universal Adapter / Barcodes (1-24)
- 1 tube PrepX PCR Primers

Package 2:

- 1 pack 384-Well Source Plate Seal (Pack of 10)

Storage Conditions

Package 1:

- Store at -20°C

Package 2:

- Store at Room temperature

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Shipping Conditions

Package 1:

- Dry ice (-70°C)

Package 2:

- Room temperature

Product Documents

Documents for our products are available for download at takarabio.com/manuals

Quality Control Data

All primers were verified by ESI analysis to have a major peak with a measured molecular weight within 0.1% of the predicted molecular weight.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: techUS@takarabio.com

United States/Canada
800.662.2566
(102518)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999

PrepX™ ILMN Barcodes 1-24 (400077-1)

CATALOG NO.

640103

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

10/25/2018