

Product Components List

ICELL8® cx Single-Cell System

| Catalog Nos. | Amount |
|--------------|--------|
| 640188 | Each |

Description

The ICELL8 cx Single-Cell System is an advanced, integrated automation platform that combines the power of imaging and dispensing with the isolation of single cells in 5,184-nanowell chips. The ICELL8 cx system employs a large-bore nozzle dispenser for unbiased cell size isolation. Imaging allows the differentiation of single cells from multiple cells. Additionally, with the use of fluorescence, imaging can also be used for live/dead analysis in conjunction with the use of cellular dyes. The full system also includes the ICELL8 cx Thermal Cycler as well as a built-in imaging station and integrated ICELL8 cx CellSelect® Software that analyzes the images and provides the ability to proceed only with chosen wells, eliminating the noise from multiple-cell-containing wells or empty wells. This open-platform, high-throughput single-cell automation system is prevalidated for use with several NGS applications and also provides the flexibility to enable users to develop applications of their choice.

Package Contents:

- 1 x ICELL8 cx Crate Components (Cat. No. 640190) [>>View Components](#)
- 5 x ICELL8 cx Loading Kit (Cat. No. 640197) [>>View Components](#)
- 5 x ICELL8 Collection Kit – L (Cat. No. 640212) [>>View Components](#)
- 1 x ICELL8 384-Well Source Plate and Seal (Cat. No. 640192) [>>View Components](#)
- 3 x ICELL8 cx Blank Chip - 250 nl FB (Cat. No. 640193) [>>View Components](#)
- 3 x ICELL8 Blank Chip Reagent Kit (Cat. No. 640196) [>>View Components](#)
- 2 x ICELL8 3' DE Reagent Kit (Cat. No. 640167) [>>View Components](#)
- 2 x ICELL8 cx 3' DE Chip (Cat. No. 640199) [>>View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- ICELL8 cx Single-Cell System User Manual
- ICELL8 cx CellSelect Software User Manual

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

| | | | |
|--------------------------------------------------|---------------------------------|------------------------------|-----------------------------|
| United States/Canada 800.662.2566 (032519) | Asia Pacific +1.650.919.7300 | Europe +33.(0)1.3904.6880 | Japan +81.(0)77.565.6999 |
|--------------------------------------------------|---------------------------------|------------------------------|-----------------------------|