

# Product Components List

---

## Shasta™ Whole-Genome Amplification Kit - 2 Chip

**Catalog No.**

640286

**Amount**

2 chips

### Description

Shasta Whole-Genome Amplification Kit - 2 Chip contains the reagents and consumables to analyze single-cell DNA sequences with whole-genome coverage using the Shasta Single Cell System (Cat. No. 640282) or ICELL8® cx Single-Cell System (Cat. No. 640188). Enough material is included for runs on two nanowell chips (5,184 wells each) following the protocol documented in the user manual.

### Package Contents

- 1 x Shasta Whole-Genome Amplification Reagents (Cat. No. 640287) (Not sold separately) [>> View Components](#)
- 2 x Shasta Long Indexing Primer Set - A (Cat. No. 640283) [>> View Components](#)
- 2 x Single-Cell 350v Chip (Cat. No. 640019) [>> View Components](#)
- 2 x Single-Cell Loading Kit (Cat. No. 640206) [>> View Components](#)
- 2 x Single-Cell Collection Kit (Cat. No. 640212) [>> View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

### Product Documents

Documents for our products are available for download at [takarabio.com/manuals](https://takarabio.com/manuals)

The following documents apply to this product:

- Shasta Whole-Genome Amplification Kit - 2 Chip User Manual for the Shasta Single Cell System
- Shasta Whole-Genome Amplification Kit - 2 Chip User Manual for the ICELL8 cx Single-Cell System

---

### Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at [takarabio.com](https://takarabio.com). It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2025 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at [takarabio.com](https://takarabio.com).

This document has been reviewed and approved by the Quality Department.

---

**Takara Bio USA, Inc.**

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

United States/Canada  
800.662.2566  
(043025)

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.565.6999