

# Product Components List

## Shasta™ Total RNA-Seq Kit - 2 Chip

Catalog No.	Amount
640288	2 chips

### Description

The Shasta Total RNA-Seq Kit - 2 Chip contains the chips, reagents, and consumables required to run two single-cell, full-length total RNA-seq experiments for up to 100,000 cells per run using the ICELL8® cx Single-Cell System (Cat. No. 640188 or 640189).

### Package Contents

- 1 x Shasta Total RNA-Seq Reagents (Cat. No. 640289) (Not sold separately) >> [View Components](#)
- 2 x Shasta Long Indexing Primer Set - A (Cat. No. 640283) >> [View Components](#)
- 2 x Single-Cell Loading Kit (Cat. No. 640206) >> [View Components](#)
- 2 x Single-Cell Collection Kit (Cat. No. 640212) >> [View Components](#)
- 2 x Single-Cell 350v Chip (Cat. No. 640019) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

### Product Documents

Documents for our products are available for download at [takarabio.com/manuals](https://takarabio.com/manuals)

The following documents apply to this product:

- Shasta Total RNA-Seq Kit - 2 Chip User Manual for the ICELL8 cx Single-Cell System

---

### Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at [takarabio.com](https://takarabio.com). It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at [takarabio.com](https://takarabio.com).

This document has been reviewed and approved by the Quality Department.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (061824)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999