

Product Components List

ThruPLEX® DNA-Seq Kit

Catalog Nos.	Amount
R400677	480 rxns

Description

The ThruPLEX DNA-Seq Kit builds on the innovative ThruPLEX chemistry to generate DNA libraries from as little as 50 pg of DNA and increases sensitivity while providing up to 96 indexes for multiplexing. This kit can be used with fragmented double-stranded DNA from any sample source—genomic DNA, biofluids such as cell-free DNA, DNA from FFPE materials, and cDNA. The entire ThruPLEX DNA-Seq Kit workflow is performed in a single tube or well in about 2 hours and requires no purification steps or sample transfers. The ThruPLEX DNA-Seq Kit can be used in any DNA-Seq application, RNA-Seq, or ChIP-Seq and offers robust target enrichment performance with all of the leading platforms.

Package Contents

- 5 x ThruPLEX DNA-Seq Kit (96 rxns) (Cat. No. R400676) [>>View component](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- ThruPLEX DNA-Seq Kit User Manual

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: technical_support@takarabio.com

United States/Canada 800.662.2566 (112219)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------