

Product Components List

Cellartis® hiPS Beta Cells (from ChiPSC22) Kit

Catalog No.	Amount
Y10106	1 kit

Description

The Cellartis hiPS Beta Cells (from ChiPSC22) Kit provides beta cells derived from a human induced pluripotent stem (hiPS) cell line. Cellartis hiPS Beta Cells (from ChiPSC22) Kit contains cryopreserved beta cells, a supplement, a coating substrate, and two types of basal media.

Package Contents

- Cellartis hiPS Beta Cells (from ChiPSC22), 1 vial (Cat. No. Y10107; not sold separately) >>[View component](#)
- Cellartis Beta Cell Supplement, 11 x 360 µl (Cat. No. Y10102; not sold separately) >>[View component](#)
- Cellartis Beta Cell Coating, 4 ml (Cat. No. Y10103; not sold separately) >>[View component](#)
- Cellartis Beta Cell Basal Medium, 160 ml (Cat. No. Y10104; not sold separately) >>[View component](#)
- Cellartis Beta Cell Basal Medium 2, 30 ml (Cat. No. Y10105; not sold separately) >>[View component](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for Takara Bio Europe AB products are available for download at www.takarabio.com/manuals

The following documents apply to this product:

- Cellartis hiPS Beta Cells Kits User Manual

Notice to Purchaser

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

This product may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Europe AB.

If you require licenses for other use, please contact us by phone at +46 31 758 0900.

Your use of this product is also subject to compliance with any applicable licensing requirements as detailed in our catalogues, on our website at <http://www.takarabio.com>, on the label or other documentation accompanying the goods. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions

Takara Bio Europe AB

A Takara Bio Company
Arvid Wallgrens backe 20, SE-413 46 Göteborg, Sweden
Europe Technical Support: techEU@takarabio.com

United States/Canada 800.662.2566 072419)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
---	---------------------------------	------------------------------	-----------------------------