

Cellartis® DEF-CS™ 500 Xeno-Free Culture Medium w/o antibiotics

Catalogue Number

Y30045

Size

1 kit

Description

Cellartis® DEF-CS™ 500 Xeno-Free Culture Medium w/o antibiotics is a pre-clinical grade, chemically defined culture medium that is free from human- and animal-derived components and is used for efficient expansion of undifferentiated human pluripotent stem cells. This kit includes basal medium and additives.

Package Contents

- Cellartis® DEF-CS™ 500 Xeno-Free Basal Medium w/o antibiotics (Y30046; not sold separately) >>[View components](#)
- Cellartis® DEF-CS™ 500 Xeno-Free Additives (Y30042) >>[View components](#)

For storage conditions, please refer to the Certificate of Analysis supplied with each component.

Product Documents

Documents for Takara Bio Europe SAS products are available for download at www.takarabio.com/manuals

The following document applies to this product:

- Cellartis® DEF-CS™ 500 Xeno-Free Culture Medium w/o antibiotics User Manual

Notice to Purchaser

This product is for Research Use Only. It is not intended for human or animal *in vivo* applications. Also, do not use this product as food, cosmetic or household item, etc.

This product may not be transferred to third parties, resold, modified for resale or transfer, or used to manufacture commercial products or to provide a service to third parties without prior written approval from Takara Bio Europe SAS.

If you require licenses for other use, please contact us by phone at +33 139 046 880.

Your use of this product is also subject to compliance with any applicable licensing requirements as detailed below, in our catalogues, on our website at <http://www.takarabio.com>, on the label or other documentation accompanying the goods. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2025 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio Europe SAS

Takara Bio Europe Swedish Filial (a Swedish branch)
Arvid Wallgrens Backe 20, SE-413 46 Göteborg, Sweden
Technical Support: techEU@takarabio.com

United States/Canada
800.662.2566
(071125)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999