

## Cellartis® iPS Cell to Hepatocyte Differentiation System

**Catalogue Number**  
Y30055

**Size**  
1 kit

### Description

The Cellartis® iPS Cell to Hepatocyte Differentiation System is a complete system for differentiation of human pluripotent stem (hPS) cells to hepatocytes via definitive endoderm (DE). This kit includes the Cellartis® DEF-CS™ Culture System (Y30020) for expansion of undifferentiated hPS cells, the Cellartis® Definitive Endoderm Differentiation Kit (Y30030) for differentiation to DE cells, and the Cellartis® Hepatocyte Differentiation Kit (Y30050) for subsequent differentiation to hepatocytes.

### Package Contents

- Cellartis® DEF-CS™ 100 Culture System (Y30020; not sold separately) >>[View components](#)
- Cellartis® Definitive Endoderm Differentiation Kit (Y30030; not sold separately) >>[View components](#)
- Cellartis® Hepatocyte Differentiation Kit (Y30050) >>[View components](#)

For storage conditions, please refer to the Certificate of Analysis supplied with each component.

### Product Documents

Documents for Takara Bio Europe SAS products are available for download at [www.takarabio.com/manuals](http://www.takarabio.com/manuals)  
The following document applies to this product:

- Cellartis® iPS Cell to Hepatocyte Differentiation System User Manual

### Notice to Purchaser

This product is for Research Use Only. It is not intended for human or animal *in vivo* applications. Also, do not use this product as food, cosmetic or household item, etc.

This product may not be transferred to third parties, resold, modified for resale or transfer, or used to manufacture commercial products or to provide a service to third parties without prior written approval from Takara Bio Europe SAS.

If you require licenses for other use, please contact us by phone at +33 139 046 880.

Your use of this product is also subject to compliance with any applicable licensing requirements as detailed below, in our catalogues, on our website at <http://www.takarabio.com>, on the label or other documentation accompanying the goods. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2025 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

---

#### Takara Bio Europe SAS

Takara Bio Europe Swedish Filial (a Swedish branch)  
Arvid Wallgrens Backe 20, SE-413 46 Göteborg, Sweden  
Technical Support: [techEU@takarabio.com](mailto:techEU@takarabio.com)

United States/Canada  
800.662.2566  
(071125)

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.565.6999

Page 1 of 1