

pDD-ZsGreen1 Reporter

Catalog No(s).

632195 (Not sold separately).
Sold as part of 632192.

Amount

20 µg pDD-ZsGreen1 Reporter Vector (500 ng/µl)

Lot Number

Specified on product label.

pDD-ZsGreen1 Reporter Vector Information




Figure 1. pDD-ZsGreen1 Reporter vector map.




Figure 2. pDD-ZsGreen1 Reporter multiple cloning site.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566
Asia Pacific +1.650.919.7300
Europe +33.(0)1.3904.6880
Japan +81.(0)77.543.6116
Late update: July 12, 2011

Location of Features

- MCS (multiple cloning site): 31–89
- DD-ZsGreen1
Start codon (ATG): 97–99; Stop codon: 1114–1116
DD (FKBP-L106P destabilization domain; 1): 97–420
ZsGreen1 (*Zoanthus* sp. green fluorescent protein): 427–1119
- f1 origin of replication (for packaging the noncoding strand of DD-ZsGreen1): 1372–1827 (complementary)
- SV40 origin of replication: 2168–2306
- Kan^r/Neo^r (kanamycin/neomycin resistance gene)
Neomycin phosphotransferase coding sequences: 2352–3146
- pUC origin of replication: 3731–4374

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.